

like never before.



2018

**AR INSTRUMED**  
DEUTSCHLAND GMBH



AR Instrumed is located in Berlin Germany having more than 2000 different product in warehouse stock available for quick delivery. The company has been retailing dental instruments and exporting them to the quality conscious healthcare industry all over the world. The company prides in its network of Germany-wide retailing and its global image of Quality instruments while dealing with international clients which repeatedly stand recognized by government as well as International business standards and Certifications, the company foresees a smiling future in dental industry.

This commitment to excellence has allowed us to establish our customer bases in strategic and vital surgical Hubs of the world, like Germany, Spain, France, Algeria.

AR INSTRUMED Deutschland GmbH has other international Business branches in its network and many distributors worldwide.

For Exclusive Distribution please contact us




## AR INSTRUMED





## Exhibition Worldwide



## Exhibition 2018-2019

IDS 2019, 13 to 17 March 2019  
Köln Germany

Many more..... Please ask

## **Supply and Delivery**

The entire range of our products is always kept available on stock in our warehouse Berlin Germany. We can thus provide you in Germany, Spain & all European with a rapid and reliable dispatch service, without middle men, directly to your practice. In export countries we supply via specialist dealers.

All orders shipped via standard shipping should be received within 5 to 7 business days. All in-stock instruments are processed and shipped within 48 hrs from our warehouse located in Berlin Germany directly after your order. You can also expedite your shipping for an extra charge by selecting priority shipping. Priority shipping orders will be arrival 1 to 2 business day after the shipping date.

## **WARRANTY POLICY**

All AR INSTRUMED products are warranted to be free from defects in workmanship and materials. Any product which proves defective in workmanship or materials will be repaired or replaced at no charge.

Additional shipping charges for use of the AR INSTRUMED DHL/Courier return labels and shipment back to the customer will be applied as shown on your invoice.

Warranty does not apply to normal wear and tear and use of the instrumentation beyond their limitations. It also excludes the improper usage of an instrument in a way that the particular instrument was designed and intended to be used. Moreover, instruments damaged in an abusive or accidental manner are also excluded by this warranty.

Warranty is also voided if the instruments are not stored and sterilized properly as shown in the Instrument Care sections of the AR INSTRUMED website and catalog. Any modification, alteration or sharpening made to AR INSTRUMED instruments will render the

## **CLEANING**

Ultrasonic cleaning is the most effective way to clean your surgical instruments, particularly instruments with hinges, box locks and other moving parts. All instruments must be pre washed prior to being fully submerged in an open position using distilled (demineralized) water and an ultrasonic cleaning solution. Always rinse ultrasonic cleaning solution from instruments with distilled water prior to sterilization.

An additional way to clean besides ultrasonic would be to manually clean the instruments with Nylon or stainless steel brushes using a mild or low PH surgical detergent. We recommend strongly to not use the following steel wool, steel brushes or scouring pads as these may damage the instrument finish. make sure to always rinse any cleaning solution from the instruments with distilled water prior to sterilization.

The common causes for staining and our suggestions on how to avoid staining of the instruments are as follows: Staining is a surface deposit on instruments, and most often mistaken for rust. After autoclaving, you may notice a stain on your instruments. Rusting instruments are very rare. Stains on instruments appear in many colors and, in most cases, the colors tell you about the origin of the stain.

**Orange/brown stain:** The problem is most often a phosphate layer (brown to light orange) on the instrument, which develops as a result of any of the following causes: water sources, detergents used to wash and clean instruments, surgical wrappings, cold sterilization solutions, or dried blood.  
**Black stain:** Black stains are commonly due to an acid reaction.

Black stains may result from detergents used to clean the instrument; similar to brown stains caused by high PH in detergents. The black acid type stain can be caused by low PH (less than six) during autoclaving.

**Dark brown stain:** Dark brown stains are usually a result of dried blood left on an instrument. Blood should be removed from the surface of the instrument immediately. It will break down the instruments surface with

## **INSTRUCTIONS FOR USE**

Dear Customer,  
We kindly ask you to read carefully the following instruction before handling and using the instruments. Each AR INSTRUMED instruments has been realized with great care and professionalism and its use is addressed only to qualified and dental specialized personnel.

### **1. FIRST USE (INITIAL USE)**

The Instrument is supplied NON-STERILED: it must properly washed, rinsed and sterilized before to be used. Many joint instruments presents oily residuals in the parts, due to the manufacturing phases; please make sure to clean the instruments into the proper way. If the use of a washing-machine is not possible, the instruments can also be washed manually, use a rigid brush with plastic bristles (No metal bristles) and make always sure the proper washing and decontamination solutions.

### **2. STERILIZATION**

It is highly recommended to sterilize the instruments only after a proper cleaning.

If the instrument is not perfectly cleaned, the result of the sterilization is prejudiced. Follow carefully the instructions for use provided by the manufacturer of the sterilization equipments.

There are different ways for sterilizing the instruments:

#### **2.1 AUTOCLAVE STERILIZATION**

The main aim of the autoclave sterilization is to obtain a proper sterilization and, in the meantime, to maintain unchanged the structural and functional characteristics of the instruments and much as possible. This process must be done only and exclusively by expert personnel, that must have been properly trained for this purpose.

A. All the cleaning operations must be carried on only in areas duly equipped and organized for this kind of operation.

B. Immediately after the use, instruments must be disinfected through the use of solutions whose effectiveness has been proved. Make sure of the right dosage of the cleaning/decontamination solutions, of the washing duration and of its temperature, following always the instructions given by the equipment manufacturer and by the cleaning solutions' producer.

C. Cleaning / disinfectant solutions containing sodium chloride (NaCl) can often cause oxidation or corrosion to the instruments compromising seriously their use.

- D. Cleaning/disinfecting solutions must be changed daily.
- E. The best recommended method for a proper washing of the instruments is the ultrasound cleaning, especially suggested for the joint instruments (ex: pliers, scissors, needle-holders,...).
- F. Do not put together in the same trays for the washing/sterilization instruments made of different metals (ex: common steel with stainless steel, aluminum, chromed steel,...).
- G. Joint instruments (ex: pliers, scissors, needle-holders,...) Must be placed in an open position for the washing disinfection.
- H. Make sure that trays are not overloaded with instruments because this could compromise the washing/disinfection.
- I. Possible residuals after the first wash must be completely removed.
- J. It has been proved that the ideal temperature for the final rinsing is between 80 C and 90 C.
- K. It is highly recommended to use only demineralized water during the final rinsing. This device avoid the rise up of dark stains.
- L. To maintain the perfect efficiency of joint instruments, it is necessary to lubricate them regularly with specific products (DO NOT USE INDUSTRIAL OIL SAT ALL).

Autoclave sterilization through saturated steam is the most used method nowadays: it is based on heat, that is necessary to sterilize, on saturated steam, that act on the thermal exchange with the instruments temperature, and on the pressure, that increase the water the water evaporation latent heat.

- For this method of sterilization, it is highly recommended to make cycles that not overcome 7 minutes, with a temperature of 134° and a pressure of 2.1 bar.
- For the enveloped instruments it is very important to complete carefully the rinsing cycle because if some water remains into the envelop, it can cause corrosion or oxidation to the instruments and generate bacteria formation.
- The sterile products must be preserved in defined areas, dry and without dust.

#### **ATTENTION: DRY HEAT STERILIZATION**

A sterilization that overcomes the 177 C (350 F) weak the temper of the instrument and consequently the cut performance and the elasticity (above all for the cut instruments, that can suffer a decrease of the sharpening); it can also occur a darkening of the instruments, brown spots that seems to oxidation.

Aluminum and Titanium: do not put aluminum instruments and titanium instruments into thermal equipment and/or ultrasound bath. For the cleaning of these instruments use only substances that are specially suggested for aluminum and titanium.

**AR INTRUMED ALL RESPONSIBILITY FOR THE FOLLOWING CASES:**

- 1. Damages caused by a not proper use of the instrument.
- 2. Damages caused by a wrong process of cleaning and/or sterilization of the instruments;
- 3. Damages caused by the warehousing of the instruments in unsuitable locals;
- 4. Damages caused by the lack of maintaining of the instruments,
- especially for the joint ones, that need a regular lubrication in the joint parts.

For further information or needs, we remain at your complete disposal. Kindly address your requests to the following address:

## **Certification**

**ISO 9001  
2008  
CERTIFIED**

**ISO 13485  
2003  
CERTIFIED**



## MANUFACTURING PROCESS

- Skilled technicians are manufacturing the shapes of each tip.
- Cutting and shaping each tip using CNC machinery.
- Proper bending of each tip, from 1 to 6 times, to insure proper use qualities.
- Increasing the strength of the instruments with the use of heat treating processes.
- AR INSTRUMED raises the level of heat treating to a 50-55 hardness rating, which increases the strength of the material and minimizes wear.
- Assembling the instruments properly is of the utmost importance. If the tip is not inserted into the handle properly, it will not work correctly and may break or fall out of the handle. Each instruments' tips are inserted by highly trained and skilled technicians.
- As a final step to the process, we clean all of the instruments and do what we call a final "blading", or sharpening, to ensure the proper working edges are on each and every instrument to meet the specifications of dental professionals.

## Handle Manufacturing Process

- Proper Cutting: To cut each handle of different sizes with precision machinery.
- Complete handle making process: To cut the tip to fit into each handle size.  
This process requires very precise work since there is a great possibility of product loss.
- Grooving work: To groove the handles so that it gives dentists comfort in treatment.
- (Grooving work is different depending on the functions and kinds of handles)
- Upscale Satin Handles: By applying friction on the handles in the sand barrel for a long time, the handles turn into bright satin finish.
- Ultrasonic Cleaning: Handles could get rusty without the ultrasonic cleaning process.
  - Without ultrasonic cleaning, it will shorten the life-span of the instrument

## ASTM STANDARD

- AR INSTRUMED applies international standards during its manufacturing. Following is a list and reference of those standards for our customers.
- Acquiring quality stainless steel as per ASTM standard F899-95 & ISO 713501.
- Corrison of Surgical Instruments in conformity to ASTM standard F1089-02
- Passivation in Conformity to ASTM standard B967-05.
- Electropolishing as per ASTM standard B.912-02.
- Rockwell Hardness test in accordance with ASTM standard E18-00
- Cleaning and Descaling of stainless steel parts as per ASTM standard A380-99 (R05)
- To Select the proper grade from AISI 410, 420 AND AISI 304 Ac, several process are undertaken according to the guidelines of ASTM standard

## Material Used in Different products

- Scissors - AISI 420 Japanese Stainless Steel.
- Gracey Curette - AISI 420 French Stainless Steel.
- Scalers - AISI 420 French Stainless Steel.
- Bone Rongeurs - AISI 420 French Stainless Steel.
- Root Elevator - TIPS AISI 420 French Stainless Steel & Handle AISI 304 Japnanese Stainless Steel..
- Chisel & Bone Scrapers - AISI 420 French Stainless Steel.
- Extracting Forceps - AISI 410 Japanese Stainless Steel.
- Forceps - AISI 410 Japanese Stainless Steel.
- Tweezers - AISI 410 Japanese Stainless Steel.
- Retractors - AISI 410 Japanese Stainless Steel.
- Needle Holder - AISI 410 Japnanese Stainless Steel.
- Orthodontics Pliers - AISI 420 Japnanese Stainless Steel.
- Rubberdam Instruments - AISI 410 Japnanese Stainless Steel.
- Mirror Handle - AISI 304 Japnanese Stainless Steel.
- Osteotomes - AISI 304 Japnanese Stainless Steel.
- Impression Tray - AISI 304 Japnanese Stainless Steel.



# Table of Content

• <i>Dental Set</i>	<i>Zahnärztliche Set</i>	<i>Conjunto dentales</i>	<i>Ensembles dentaires</i>	008 to 024
• <i>Extraction</i>	<i>Extraktion</i>	<i>Extracción</i>	<i>Extraction</i>	025 to 063
• <i>Diagnostic</i>	<i>Diagnostik</i>	<i>Diagnóstico</i>	<i>Diagnostique</i>	064 to 069
• <i>Periodontology</i>	<i>Parodontologie</i>	<i>Periodoncia</i>	<i>Parodontologie</i>	070 to 078
• <i>Conservative</i>	<i>Konservierung</i>	<i>Conservador</i>	<i>Soins Conservateur</i>	079 to 088
• <i>Surgery</i>	<i>Chirurgie</i>	<i>Cirugía</i>	<i>Chirurgie</i>	089 to 121
• <i>Micro Surgery</i>	<i>MikroChirurgie</i>	<i>Micro Cirugía</i>	<i>Microchirurgie</i>	122 to 125
• <i>Impression Trays</i>	<i>Abdrucklöffel</i>	<i>Cubetas de Impresión</i>	<i>Portes empreintes</i>	126 to 130
• <i>Implantology</i>	<i>Implantologie</i>	<i>Implantología</i>	<i>Implantologie</i>	131 to 144
• <i>Orthodontics</i>	<i>Kieferorthopädie</i>	<i>Ortodoncia</i>	<i>Orthodontie</i>	145 to 167
• <i>Prosthesis</i>	<i>Prothese</i>	<i>Prótesis</i>	<i>Prothèse</i>	168 to 173
• <i>Calipers</i>	<i>Messinstrumente</i>	<i>Calibres</i>	<i>Instruments de mesure</i>	174 to 176
• <i>Laboratory</i>	<i>Laboratorium</i>	<i>Laboratorio</i>	<i>Laboratoire</i>	177 to 180
• <i>Sterilization</i>	<i>Sterilisation</i>	<i>Esterilización</i>	<i>Stérilisation</i>	181 to 185
• <i>Dental Implant</i>	<i>Zahnimplantat</i>	<i>Implante dental</i>	<i>Implant dentaire</i>	186 to 195



*Dental Set - Zahnärztliche Set - Conjunto dentales - Ensembles dentaires*



**AR INSTRUMED**

Extracting Forceps Sets  
Zahnzangen Sets  
Fórceps para extracción Juegos  
Davier dentaire pour extraction Set



**110-006**  
Extracting forceps  
Set for adults  
set of 10 pcs.



**110-048**  
Micro  
Extracting forceps  
Set for adults  
set of 10 pcs.



**110-007**  
Extracting forceps  
Set for childrens  
set of 7 pcs.

Fig 1:100-001 Upper incisors & Canines  
Fig 7:100-005 Upper Premolars  
Fig 13:100-007 Lower Premolars  
Fig 17:100-010 Upper Molars, Right  
Fig 18:100-011 Upper Molars, Left  
Fig 22:100-017 Lower Molars  
Fig 33A:100-029 Lower Roots  
Fig 51:100-042 Upper Roots  
Fig 67:100-054 Upper Third Molars  
Fig 79A:100-071 Lower Wisdoms

Klein 100-101 Upper Incisors  
Klein 100-102 Upper Premolars  
Klein 100-103 Upper Molars  
Klein 100-104 Upper Roots  
Klein 100-105 Lower Incisors  
Klein 100-106 Lower Molars  
Klein 100-107 Lower Roots



**AR INSTRUMED**



Extracting Forceps Sets  
 Zahnzangen Sets  
 Fórceps para extracción Juegos  
 Davier dentaire pour extraction Set

Fig 1:100-001A Upper incisors & Canines  
 Fig 7:100-005A Upper Premolars  
 Fig 13:100-007A Lower Premolars  
 Fig 17:100-010A Upper Molars, Right  
 Fig 18:100-011A Upper Molars, Left  
 Fig 22:100-017A Lower Molars  
 Fig 33A:100-029A Lower Roots  
 Fig 51:100-042A Upper Roots  
 Fig 67:100-054A Upper Third Molars  
 Fig 79A:100-071A Lower Wisdoms

110-018 Black Plasma  
 Extracting forceps  
 Set for adults  
 set of 10 pcs  
 Special Black  
 Plasma coating

110-031A  
 Extracting forceps  
 Set for adults  
 set of 11 pcs.



**AR INSTRUMED**

Extraction & Luxators Sets  
Extraktion & Wurzelheber Sets  
Extracción & Luxadores Sets  
Extraction & Élèvateurs Sets



104-062

Root Elevator Luxators  
Set of 7

- 104-063 Fig. 1 Root Elevator luxatur 2 mm
- 104-064 Fig. 2 Root Elevator luxatur 3 mm
- 104-065 Fig. 3 Root Elevator luxatur 5 mm
- 104-066 Fig. 4 Root Elevator luxatur 3 mm
- 104-067 Fig. 5 Root Elevator luxatur 5 mm
- 104-068 Fig. 6 Root Elevator luxatur 3 mm
- 104-069 Fig. 7 Root Elevator luxatur 3 mm



110-121

Extraction Elevator  
Set of 10



**AR INSTRUMED**

Extraction & Luxators Sets  
Extraktion & Wurzelheber Sets  
Extracción & Luxadores Sets  
Extraction & Eléveurs Sets



101-301

- |                |                                      |
|----------------|--------------------------------------|
| Luxators       | 101-301-1 Luxator 1mm Straight       |
| Plastic Handle | 101-301-2 Luxator 2mm Straight       |
| Set of 13      | 101-301-3 Luxator 3mm Straight       |
|                | 101-301-4 Luxator 3mm Curved         |
|                | 101-301-5 Luxator 3mm Curved Upward  |
|                | 101-301-6 Luxator 2mm Curved         |
|                | 101-301-7 Luxator 4mm Curved         |
|                | 101-301-8 Luxator 5mm Straight       |
|                | 101-301-9 Luxator 5mm Curved         |
|                | 101-301-10 Luxator 3mm/5mm           |
|                | 101-301-11 Luxator 2mm/3mm           |
|                | 101-301-12 Luxator 3mm Double Curved |
|                | 101-301-13 Luxator 3mm Straight      |



110-034  
Luxators  
Anatomical Handle  
Set of 7



110-047  
Heidbrink  
Set of 3  
Round Handle

- 101-128 Fig. 1 Heidbrink Straight 1.8mm  
101-128 Fig. 2 Heidbrink Right 1.8mm  
101-128 Fig. 3 Heidbrink Left 1.8mm



**AR INSTRUMED**

Extraction & Micro Surgery  
Extraktion & Mikro Chirurgie  
Extracción & Micro Cirugía  
Extraction & Micro Chirurgie



104-052A

Extraction Instruments  
Set of 9

- 100-001 Fig. 1 Extracting Forcep
- 100-042 Fig. 51 Extracting Forcep
- 100-017 Fig. 22 Extracting Forcep
- 100-005 Fig. 7 Extracting Forcep
- 101-097 Fig. 2 Elevator Bein 3 mm
- 101-014 Elevator Flohr Right 3 mm
- 101-015 Elevator Flohr Left 3 mm
- 103-140 Grip Needle Holder Mathieu 14 cm TC
- 103-131-1 Grip Needle Holder Mayo Hegar 15 cm TC



109-100

Micro Surgery  
Set of 10

- 109-006 Micro Iris Scissors Curved 105mm Blade 18mm TC Serration
- 109-007 Micro Iris Scissors Straight 105mm Blade 30mm TC Serration
- 109-008 Micro Iris Scissors Curved 105mm Blade 30mm TC Serration
- 109-009 Micro Iris Scissors Curved 105mm Blade 15mm TC Serration
- 109-010 Micro Iris Scissors Curved 130mm Blade 15mm TC Serration
- 109-011 Micro Iris Scissors Curved 105mm Blade 20mm TC Serration
- 109-012 Micro Iris Scissors Curved 105mm Blade Round TC
- K25-5006 Micro Castroviejo Neede Holder Curved 140mm TC
- 103-076 Micro Castroviejo Scissors Curved 140mm
- 102-099 TC Round Tweezer 140mm TC



**AR INSTRUMED**



- 106-339 Bone Scraper
- 106-340 Chisel 1 mm / 4 mm
- 106-341 Chisel 1.5 mm / 4 mm
- 106-342 Chisel 2 mm / 4 mm
- 106-343 Chisel 2.5 mm / 4 mm
- 106-362 Bone Chisel 5 mm Centre
- 106-366 Bone Chisel 5 mm Left
- 106-364 Bone Chisel 5 mm Right
- 106-059H Fig. 1 Sinus Lift Instrument
- 106-059H Fig. 2 Sinus Lift Instrument
- 106-059H Fig. 3 Sinus Lift Instrument
- 106-059H Fig. 4 Sinus Lift Instrument

110-021

Implantology  
Set of 12



- 104-401 William Periodontal Probe
- 104-402 Goldman-Fox Universal Scaler Fig. 1
- 104-403 Mc Call Scaler Fig. 13/14
- 104-404 Crane-Kaplan Scaler Fig. 6
- 104-405 Towner Scaler Fig. 15
- 104-406 Gracey Curette Fig.5/6
- 104-407 Gracey Curette Fig.11/12
- 104-408 Gracey Curette Fig.13/14**
- 104-409 Kirkland Gingivectomy Knife Fig. 15k/16k
- 104-410 Orban Gingivectomy Knife Fig.1/2
- 104-411 Buck Gingivectomy Knife Fig. 5/6
- 104-412 Molt Periosteal Elevator Fig. M9
- 104-413 Goldman Fox Scissors Curved
- 104-414 Crane Kaplan Right
- 104-415 Crane Kaplan Left



**AR INSTRUMED**

Rubber dam Clamps

Kofferdam

Porta Clamps

Set de Pince Porte Clamp/Pinces à Digue



110-025

Rubber dam Punches  
Set of 4

103-271 Rubber Dam Punch ANISWORTH

103-272 Rubber Dam Punch IVORY

103-274 Rubber Dam Punch IVORY LIGHT WEIGHT

103-276 Rubber Dam Punch STOCKES



110-025-1

Rubberdam Kit  
17 Pcs  
Leather Kit




110-025-2

Rubberdam Kit  
12 Pcs  
Plastic Box



**AR INSTRUMED**



Implantology Set  
Implantology Set  
Implantology Set  
Implantology Set

102-004	Mirror Plane 4
102-015-1	Mirror Handle Hollow 10mm
100-322T	Suction tube Titanium
100-241	Friedman 3.5mm 14cm
102-032	Adson Brown 7:7 Teeth
102-067	Mariam Tweezer
103-050-1	Iris Scissors 11.5cm Curved TC
103-090	Lagrange Scissors 14cm
103-134-1	Needle Holder Crile Wood 15cm TC
103-165	Halstead Mosquito Curved 12cm
PI-336	Towel Forceps with Clip
106-217	Castroviejo Caliper 20mm
106-247	Minnesota Retractor
MH-366	Scalpel Handle Scale
MH-367	Scalpel Handle Round Straight
106-131	Seldin 23
102-902H	Gracey Curette 13-14
106-051H	Prichard Periosteal Elevator
103-420	Molt 9 Periosteal Elevator
103-466	Surgical Curette 2.5mm
103-469	Surgical Curette 4mm
106-050	Periosteal Elevator P24G
103-415	Mead Hammer 17cm
106-362G	Bone Chisel Straight 5mm
106-275	Basins for Bone

110-032

Implant Instruments  
Set of 25  
Stainless Steel Tray



**AR INSTRUMED**

Perforated  
Perforierte  
Perforada  
Solido  
Perfore  
Solid  
Solide  
Solida  
Plein



Set of 6  
107-350  
S-M-L  
Upper/Lower

**Feature**  
No Welding  
Easy to Clean  
Easy Sterilization/Autoclave  
Sterilizing Autoclave at 134c  
Adjustable Front Wall  
High Quality AISI 304 Stainless Steel  
Micro Sand Blast Finished/Mirror Finished



Set of 6  
107-370  
S-M-L  
Upper/Lower

**Feature**  
No Welding  
Easy to Clean  
Easy Sterilization/Autoclave  
Sterilizing Autoclave at 134c  
Adjustable Front Wall  
High Quality AISI 304 Stainless Steel  
Micro Sand Blast Finished/Mirror Finished



**AR INSTRUMED**

Perforated  
Perforierte  
Perforada  
Solido  
Perfore  
Solid  
Solide  
Solida  
Plein



**Set of 6**  
107-390  
Full Denture  
Perforated  
S-M-L  
Upper/Lower

**Feature**  
No Welding  
Easy to Clean  
Easy Sterilization/Autoclave  
Sterilizing Autoclave at 134c  
Adjustable Front Wall  
High Quality AISI 304 Stainless Steel  
Micro Sand Blast Finished/Mirror Finished



**Set of 6**  
107-410  
Full Denture  
Solid  
S-M-L  
Upper/Lower

**Feature**  
No Welding  
Easy to Clean  
Easy Sterilization/Autoclave  
Sterilizing Autoclave at 134c  
Adjustable Front Wall  
High Quality AISI 304 Stainless Steel  
Micro Sand Blast Finished/Mirror Finished



**AR INSTRUMED**

Solid  
Solide  
Solido  
Plein

**Set of 6**  
107-037-2  
S-M-L  
Upper/Lower



**Set of 8**  
107-137  
S-M-L-XL  
Upper/Lower



**Impression Tray  
Rimlock  
Solid**



**Set of 10**  
107-137-1  
XS-S-M-L-XL  
Upper/Lower



**AR INSTRUMED**

Perforated  
Perforierte  
Perforada  
Perfore

Impression Tray  
Rimlock  
Perforated



**Set of 6**  
107-136  
S-M-L  
Upper/Lower

**Set of 8**  
107-133  
S-M-L-XL  
Upper/Lower



**AR INSTRUMED**

Solid / Perforated  
Solide / Perforierte  
Solido / Perforada  
Plein / Perfore

Set of 6  
75-386  
S-M-L  
Upper/Lower



Impression Tray  
Rimlock  
Perforated  
Solid

Set of 8  
107-092-1  
S-M-L-XL  
Upper/Lower

Set of 8  
107-106-1  
S-M-L-XL  
Upper/Lower

Gracey Curettes



110-005  
Gracey Curettes  
Set of 7  
10mm  
Hollow Handle



**AR INSTRUMED**

Dental Set  
Zahnärztliche Set  
Conjunto dentales  
Ensembles dentaires



110-117

Examination Set of 4  
Explorer  
Mirror  
Mirror Handle  
Tweezer



110-118

Examination Set of 4  
Explorer  
Mirror  
Mirror Handle  
Tweezer



110-110

Examination Set of 4  
Explorer  
Mirror  
Mirror Handle  
Tweezer



110-118-1

Examination Set of 5  
Explorer  
Mirror  
Mirror Handle  
Tweezer  
Heidmann Spatula



75-379

Gracey Curettes  
Set of 8  
Leather Case  
6mm Handle

*Dental Set - Zahnärztliche Set - Conjunto dentales - Ensembles dentaires*



**AR INSTRUMED**

Dental Set  
Zahnärztliche Set  
Conjunto dentales  
Ensembles dentaires

75-376

PK Thomas  
Set of 5  
6mm Handle  
Leather Case

110-119

Composite Instruments  
Set of 6  
10mm Hollow Handle  
Gold Plasma  
Leather Case

104-416A

Instruments  
Set of 13  
Leather Case

110-115

Laboratory Instruments  
Set of 16  
Leather Case

*Dental Set - Zahnärztliche Set - Conjunto dentales - Ensembles dentaires*



**AR INSTRUMED**

Extracting Forceps Sets  
Zahnzangen Sets  
Fórceps para extracción Juegos  
Davier dentaire pour extraction Set

110-030A

Extracting Forceps  
Set of 11  
Leather Case



Fig 1:100-001A Upper incisors & Canines  
Fig 7:100-005A Upper Premolars  
Fig 13:100-007A Lower Premolars  
Fig 17:100-010A Upper Molars, Right  
Fig 18:100-011A Upper Molars, Left  
Fig 22:100-017A Lower Molars  
Fig 33A:100-029A Lower Roots  
Fig 51:100-042A Upper Roots  
Fig 67:100-054A Upper Third Molars  
Fig 79A:100-071A Lower Wisdoms  
Fig 86C : 100-076A Lower Molar

Fig 1:100-001 Upper incisors & Canines  
Fig 7:100-005 Upper Premolars  
Fig 13:100-007 Lower Premolars  
Fig 17:100-010 Upper Molars, Right  
Fig 18:100-011 Upper Molars, Left  
Fig 22:100-017 Lower Molars  
Fig 33A:100-029 Lower Roots  
Fig 51:100-042 Upper Roots  
Fig 67:100-054 Upper Third Molars  
Fig 79A:100-071 Lower Wisdoms  
Fig 86C : 100-076 Lower Molar



110-031  
Anatomical  
Extracting Forceps  
Set of 11  
Leather Case



*Extraction-Extraktion-Extracción-L'Extraction*



**AR INSTRUMED**

Extracting Forceps

Zahnzangen

Fórceps para extracción

Davier dentaire pour extraction





**AR INSTRUMED**

Extracting Forceps  
Zahnzangen  
Fórceps para extracción  
Davier dentaire pour extraction





**AR INSTRUMED**

Extracting Forceps  
Zahnzangen  
Fórceps para extracción  
Davier dentaire pour extraction





**AR INSTRUMED**

Extracting Forceps  
Zahnzangen  
Fórceps para extracción  
Davier dentaire pour extraction





**AR INSTRUMED**

Extracting Forceps  
Zahnzangen  
Fórceps para extracción  
Davier dentaire pour extraction





**AR INSTRUMED**

Extracting Forceps  
Zahnzangen  
Fórceps para extracción  
Davier dentaire pour extraction





**AR INSTRUMED**

Extracting Forceps  
Zahnzangen  
Fórceps para extracción  
Davier dentaire pour extraction





# AR INSTRUMED

Black Scorpio and extraction forceps for removal of the fragments, after resection and residues of the teeth eightths inferior included especially if mesioverso this may also be used for the upper jaw in the case of difficult access to the roots. The particular shape of the forcep to avoid interference with the adjacent teeth, and prevents so to leverage accidentally on the seventh & on the sixth of arch concerned.

Made of high quality steel AISI 420, the instrument becomes characteristics of corrosion and coated with a special surface treatment, the ergonomic and lightweight slender beaks me but from the power decided will give the professional a strong feeling of use pliers solid, powerful yet precise and safe.

Therefore it will be considerable and determined also the reduction of time d intervention and consequently the stress of the patient and the operator's.



Black  
Scorpion

100-097

Forcep  
for removing  
fragments high  
quality steel  
AISI 420  
170mm



100-098

Fig 1  
Meissner Forcep  
3mm

100-099

Fig 2  
Meissner Forcep  
4mm

100-100

Fig 3  
Meissner Forcep  
6.5mm



**AR INSTRUMED**

**A-TRAUMATIC EXTRACTION SET**

Improved tapered beaks for complete root contact.

Cross Serration for a never slip grasp.

Easily Enter the sub gingival.

NO FRACTURE-NO TISSUE DAMAGE-EXTRACTION WITH EASE

Extracting Forceps

Zahnzangen

Fórceps para extracción

Davier dentaire pour extraction



75-352

Upper Anteriors  
(Adults)



75-348

Upper Universal  
(Adults)



75-349

Lower Universal  
(Adults)



75-351

Lower Molars  
(Adults)



75-353

Lower Anteriors  
(Adults)



# AR INSTRUMED

**A-TRAUMATIC EXTRACTION SET**  
Improved tapered beaks for complete root contact.  
Cross Serration for a never slip grasp.  
Easily Enter the sub gingival.  
NO FRACTURE-NO TISSUE DAMAGE-EXTRACTION WITH EASE

Extracting Forceps  
Zahnzangen  
Fórceps para extracción  
Davier dentaire pour extraction





**AR INSTRUMED**

Extracting Forceps  
Zahnzangen  
Fórceps para extracción  
Davier dentaire pour extraction





**AR INSTRUMED**

Extracting Forceps  
Zahnzangen  
Fórceps para extracción  
Davier dentaire pour extraction



200-21  
lower Molars  
17cm



200-8  
lower Premolars  
16.5cm



200-130  
Forcep for  
sevenths - eighths  
superior 19.5cm




200-187  
Forcep for  
fragments of  
lower roots  
19cm





**AR INSTRUMED**


Extracting Forceps  
Zahnzangen  
Fórceps para extracción  
Davier dentaire pour extraction





**AR INSTRUMED**

Extracting Forceps  
Zahnzangen  
Fórceps para extracción  
Davier dentaire pour extraction





**AR INSTRUMED**

Extracting Forceps  
Zahnzangen  
Fórceps para extracción  
Davier dentaire pour extraction



100-058

Fig 73  
lower Molars



100-064

Fig 74  
Lower Roots



100-081

Fig 99  
Lower Molars



100-071

Fig 79 A  
lower Wisdoms



100-076

Fig 86C  
lower Molar  
Cow Horn Beak



**AR INSTRUMED**

Children Extracting Forceps  
Kinder Zahnzangen  
Fórceps para extracción niños  
Enfants Davier dentaire pour extraction



100-101  
Klein  
Upper Incisors



100-102  
Klein  
Upper Premolars



100-103  
Klein  
Upper Molars



100-104  
Klein  
Upper Roots



100-105  
Klein  
Lower Incisors



100-106  
Klein  
Lower Molars



100-107  
Klein  
Lower Roots



**AR INSTRUMED**

Extracting Forceps  
Zahnzangen  
Fórceps para extracción  
Davier dentaire pour extraction



**100-230**

Root Forceps  
45 Angled



**100-231**

Root Forceps  
90 Angled



**100-230-1**

Root Forceps  
Anatomica Handle  
45 Angled



**AR INSTRUMED**

Forceps  
Klemmen  
Pinzas  
Pinces



102-114

Forcep for  
fragments of root  
Straight  
15cm



102-115

Forcep for  
fragments of root  
Angled 45°  
15cm



**AR INSTRUMED**

Bone Rongeurs  
Hohlmeißelzangen  
Pinzas gubias  
Pinces-Gouges





**AR INSTRUMED**

Bone Rongeurs  
Hohlmeißelzangen  
Pinzas gubias  
Pinces-Gouges



100-274

Cleveland  
Bone Rongeur  
14cm

100-244

Micro-Friedman  
Bone Rongeur  
16cm

100-243

Micro-Friedman  
Bone Rongeur  
S-Shaped  
15cm

100-247

Micro-Friedman  
Bone Rongeur  
14cm  
-90°



**AR INSTRUMED**

Bone Rongeurs  
Hohlmeißelzangen  
Pinzas gubias  
Pinces-Gouges



100-238

Blumenthal  
Bone Rongeur  
15cm  
-90°

100-239

Blumenthal  
Bone Rongeur  
15cm  
45°

100-240

Blumenthal  
Bone Rongeur  
15cm  
30°

100-241

Friedman  
Bone Rongeur  
14cm



100-329-1

Intraligamentary  
syringe  
Stainless Steel  
1.8ml  
EU needle  
Blue Plasma

The anesthesia intraligamentary allows to inject the anesthetic directly into the ligament of the tooth so that ensure the rapid and total anesthesia only the affected tooth without having to anesthetize instead more generically an area as happens with normal anesthesia

### Product Features

Tip angle of 45 ° in order to make more convenient the operativity.

Insert anesthetic with finest

Handle with lever dosage convenient to work in both the top and the bottom of the oral cavity

Processing much less traumatic than to classical anesthesia.

Less quantities of anesthetic given for each dose 0.06 ml for single push



100-329

Intraligamentary  
syringe  
Stainless Steel  
1.8ml  
EU needle

### Advantages

Quick effect of anesthetic

Eliminated the risk hematoma formation

Plus possibility of processing of teeth located in different position.

Effect of anesthetic circumscribed even just at a single element to be treated.

Easy of operate.



**AR INSTRUMED**

Syringes  
Spritzen  
Jeringas  
Seringues

*Intraligamentary Syringe*



100-288  
Intraligamentary  
syringe  
1.8ml (0,07 ml/ click)  
EU needle

*Intraligamentary Syringe*



100-284  
Intra-ligamental  
syringe  
1.8ml  
EU needle



100-311  
T Handle Type  
syringe  
1.8ml  
EU needle



100-302  
Single ring type  
syringe  
1.8ml  
EU needle



**AR INSTRUMED**

Syringes  
Spritzen  
Jeringas  
Seringues


*Syringes for Anesthesia*



100-297  
Single Ring Type  
syringe  
1.8ml  
EU needle



100-313  
T- handle type  
syringe  
1.8ml  
EU needle



100-305  
3 Ring type  
syringe  
1.8ml  
EU needle



100-298  
3 Ring type  
syringe  
1.8ml  
EU needle



# AR INSTRUMED

Extraction Elevators  
Extraktion Wurzelheber  
Extracción Elevadores  
Extraction Elévateurs

The new line of tips Micro Elevators of AR INSTRUMED easily to operate because the insertion of the lever in intraligamentary will facilitated thanks to ergonomics and accuracy in all essential shape.

All this has been studied in order to ensure less traumatic during the extraction and consequently faster post-operative recoveries.

The tips of the instruments are treated with a special coating of titanium allowing the use even for implant.





**AR INSTRUMED**

Root Elevators  
Wurzelheber  
Elevadores  
Élèvateurs



101-097

Fig 2  
Bein 3 mm




101-098

Fig 3  
Bein 4 mm



101-155

Fig 3  
Bayonet distal  
4mm



101-157

Fig 2  
Bayonet Mesial  
3mm




101-021

Fig 1  
Lindo-Levien  
2.5 mm



101-022

Fig 2  
Lindo-Levien  
3.5 mm



101-023

Fig 3  
Lindo-Levien  
4.5 mm



**AR INSTRUMED**

Root Elevators  
Wurzelheber  
Elevadores  
Élèvateurs





**AR INSTRUMED**

Root Elevators  
Wurzelheber  
Elevadores  
Élèvateurs



101-128

Fig 1  
Root - Picker  
1.8 mm  
15cm



101-129

Fig 2  
Root - Picker  
1.8 mm  
15cm



101-130

Fig 2  
Root - Picker  
1.8 mm  
15cm



430 - 8

Heidbrink  
2.5 mm  
15cm



430 - 9

Heidbrink  
2.5 mm  
15cm



430 - 10

Heidbrink  
2.5 mm  
15cm



101-128-1

H 1  
Heidbrink  
1.8 mm  
Single End



101-129-1

H 2  
Heidbrink  
1.8 mm  
Single End



101-130-1

H 3  
Heidbrink  
1.8 mm  
Single End



**AR INSTRUMED**

Root Elevators  
Wurzelheber  
Elevadores  
Élèvateurs



101-204

Fig 1  
Chompret



101-206

Fig 3  
Chompret



101-207

Fig 4  
Chompret



104-198

Warwick-James  
Straight  
3mm



104-199

Warwick-James  
Left  
6.5mm



104-200

Warwick-James  
Right  
6.5mm



101-067

Fig 25  
Cryer



101-068

Fig 26  
Cryer

54

*Extraction-Extraktion-Extracción-Extraction*



**AR INSTRUMED**

Root Elevators  
Wurzelheber  
Elevadores  
Élèvateurs



480-1

Seldin 3mm  
Straight



480-2

Seldin 4.5mm  
Straight



480-3

Hylin 3mm



480-4

Hylin 3mm



480-5

Hylin special  
for molars and  
premolars  
4mm



**AR INSTRUMED**

Root Elevators  
Wurzelheber  
Elevadores  
Élèvateurs



505-1  
Kopp 3mm  
13cm



505-2  
Kopp 3mm  
13cm



505-3  
Kopp 4mm  
13cm



505-4  
Kopp 4mm  
13cm



505-5  
Kopp 4mm  
12.5cm



505-6  
Kopp 4mm  
12.5cm



**AR INSTRUMED**

Root Elevators  
Wurzelheber  
Elevadores  
Élèvateurs



101-233

Luxator  
Straight  
2 mm

104-063

Luxator  
Cvd  
Backwad  
2 mm



101-234

Luxator  
Straight  
2.5 mm



101-235

Luxator  
Straight  
3 mm

104-066

Luxator  
Cvd  
Backwad  
3 mm

104-068

Luxator  
Cvd  
Forward  
3 mm

104-069

Luxator  
Double  
Cvd  
3 mm



101-236

Luxator  
Straight  
3.5mm



101-237

Luxator  
Straight  
4 mm



101-238

Luxator  
Straight  
5 mm

104-067

Luxator  
Cvd  
Backwad  
5 mm



**AR INSTRUMED**

Root Elevators  
Wurzelheber  
Elevadores  
Élèvateurs





**AR INSTRUMED**

Root Elevators  
Wurzelheber  
Elevadores  
Élèvateurs



75-043	75-044	75-045	75-046	75-047	75-048	75-049
Elevator straight head 2 mm Finger hollow handle	Elevator straight angled head 3 mm Finger hollow handle	Elevator straight blade head 3 mm Finger hollow handle	Elevator Contra angle head 3 mm Finger hollow handle	Elevator Inverted Curved head 3 mm Finger hollow handle	Elevator Curved Blade internal 5 mm Finger hollow handle	Elevator Straight Blade external 5 mm Finger hollow handle



**AR INSTRUMED**

Root Elevators  
Wurzelheber  
Elevadores  
Élèvateurs



101-180

Walter F Barry

101-181

Walter F Barry

101-172

Root Elevator

101-173

Root Elevator

101-182

Winter F Barry

101-183

Winter F Barry



**AR INSTRUMED**

Root Elevators  
Wurzelheber  
Elevadores  
Élèvateurs

Plastic Handle Luxators / Sterilizable



101-246	101-247	101-248	101-249	101-250	101-251	101-252	101-253	101-254	101-255	101-256
Luxator 2mm Straight	Luxator 2mm Curved	Luxator 3mm Straight	Luxator 3mm Curved	Luxator 5mm Straight	Luxator 5mm Curved	Luxator 2mm Curved Backward Front side	Luxator 2mm Curved	Luxator 1.5/3mm Straight	Luxator 2/3mm Straight	Luxator 3/5mm Straight



**AR INSTRUMED**

Root Elevators  
Wurzelheber  
Elevadores  
Élèvateurs



**101-240**

Extraction Elevator  
2mm  
Right

**101-241**

Extraction Elevator  
2mm  
Left

**101-242**

Extraction Elevator  
3mm  
Right

**101-243**

Extraction Elevator  
3mm  
Left

**101-244**

Extraction Elevator  
3mm  
Right

**101-245**

Extraction Elevator  
3mm  
Left



**AR INSTRUMED**

Periotomes  
Periotom  
Periotomes  
Periotomes



EPI-2S

Periotome  
Flexibe  
2mm  
Straight

EPI-2C

Periotome  
Flexibe  
2mm  
Curved

EPI-3S

Periotome  
Flexibe  
3mm  
Straight

EPI-3C

Periotome  
Flexibe  
3mm  
Curved

EPI-3M2C

Periotome  
Flexibe  
2mm  
Curved

EPI-3MC

Periotome  
Flexibe  
2mm  
Curved



*Diagnostic-Diagnostik-Diagnóstico-Diagnostic*

---



# AR INSTRUMED


Mouth mirrors handle  
Mundspiegelgriff  
Mango de espejos  
Manche de miroir





# AR INSTRUMED

Mouth mirrors  
Spiegel  
Espejos  
miroir



102-004	102-007	102-010	102-006	102-009	102-012	102-005	102-008	102-004-1	102-007-1	102-006-1	102-009-1
Mouth mirror No. 4 plane Ø 22mm Universal	Mouth mirror No. 5 plane Ø 24mm Universal	Mouth mirror No. 6 plane Ø 26mm Universal	Mouth mirror No. 4 Front-Surface Ø 22mm Universal	Mouth mirror No. 5 Front-Surface Ø 24mm Universal	Mouth mirror No. 6 Front-Surface Ø 26mm Universal	Mouth mirror No. 4 Magnifying Ø 22mm Universal	Mouth mirror No. 5 Magnifying Ø 24mm Universal	Mouth mirror No. 4 plane Ø 22mm Cone-Socket	Mouth mirror No. 5 plane Ø 24mm Cone-Socket	Mouth mirror No. 4 Front Surface Ø 22mm Cone-Socket	Mouth mirror No. 5 Front Surface Ø 24mm Cone-Socket





**AR INSTRUMED**

Tweezers  
Pincetten  
Pinzas  
Précelle



102-059  
Meriam  
Serrated  
16cm



102-074  
London-College  
Serrated  
15cm



102-087  
Self-Locking  
Serrated  
15cm



102-067  
Meriam  
Serrated  
16cm



**AR INSTRUMED**

Explorers  
Sonden  
Exploradores  
Sondes



102-127H  
Explorer  
130 mm



102-130H  
Explorer  
130 mm



102-124H  
Explorer  
130 mm



102-120H  
Explorer



102-125  
Explorer

102-119  
Explorer

102-120  
Explorer

102-127  
Explorer

Root-Canal Explorers with anatomical hollow handle ø 10 mm

Root-Canal Explorers with anatomical hollow handle ø 6 mm



# AR INSTRUMED

Explorers  
Sonden  
Exploradores  
Sondes



Root-Canal Explorers with anatomical hollow handle ø6 mm



Spreaders for lateral condensation with anatomical hollow handle ø6 mm



*Periodontology-Parodontologie-Periodoncia-Parodontologie*

---



# AR INSTRUMED

Periodontal Probes  
Parodontal Sonden  
Sondas Parodontales  
Sondes parodontales




Hollow handle 10mm

Hollow handle 6mm



**AR INSTRUMED**

Curettes Gracey  
Gracey Küretten  
Curetas Gracey  
Curettes de gracey




Curettes Gracey with Zirconium Nitride(ZrN) coating tips with anatomical hollow handle 10 mm



**AR INSTRUMED**

Curettes Gracey  
Gracey Küretten  
Curetas Gracey  
Curettes de gracey



102-896H / 1-2  
Curette Gracey  
Anatomical hollow  
handle ø10 mm  
175mm

102-897H / 3-4  
Curette Gracey  
Anatomical hollow  
handle ø10 mm  
175mm

102-898H / 5-6  
Curette Gracey  
Anatomical hollow  
handle ø10 mm  
175mm


102-899H / 7-8  
Curette Gracey  
Anatomical hollow  
handle ø10 mm  
175mm

102-900H / 9-10  
Curette Gracey  
Anatomical hollow  
handle ø10 mm  
175mm

102-901H / 11-12  
Curette Gracey  
Anatomical hollow  
handle ø10 mm  
175mm

102-902H / 13-14  
Curette Gracey  
Anatomical hollow  
handle ø10 mm  
175mm

102-903H / 15-16  
Curette Gracey  
Anatomical hollow  
handle ø10 mm  
175mm



Curettes Gracey with anatomical hollow handle 10 mm



**AR INSTRUMED**


Curettes Gracey  
Gracey Küretten  
Curetas Gracey  
Curettes de gracey





# AR INSTRUMED

Set of implant cleaning scalers & curettes  
Set zur reinigung von implantaten  
Set de curetas y limpiadores para implantes  
Set de Curettes pour Implantologie



104-405  
Towner Scaler  
Fig. 15

Periodontal scalers kit with anodization of Titanium (Blue) coating tips with anatomical hollow handle 6 mm



**AR INSTRUMED**

Scaler  
Scaler  
Scaler  
Scaler



102-178H

Scaler Mitchell  
175 mm



102-185H

Scaler



106-127H

Scaler Orban  
Knife



106-128H

Scaler Krikland  
Knife



102-193H

Scaler



102-184H

Scaler Carne  
Kaplan



Anatomical hollow handle 10 mm



**AR INSTRUMED**

Scaler  
Scaler  
Scaler  
Scaler



102-178  
Mitchell  
osteotrimmer  
scaler  
anatomical  
hollow  
handle  
6 mm



102-191  
Towner  
Jacquette  
scaler  
anatomical  
hollow  
handle  
6 mm



102-951  
H6/H7  
scaler  
anatomical  
hollow  
handle  
6 mm



102-234  
Kirkland  
scaler  
anatomical  
hollow  
handle  
6 mm



102-192  
Towner  
scaler  
anatomical  
hollow  
handle  
6 mm



102-184  
Crane kaplan  
scaler  
anatomical  
hollow  
handle  
6 mm



102-451  
Hollen back  
1/2scaler  
anatomical  
hollow  
handle  
6 mm



**AR INSTRUMED**

Scaler  
Scaler  
Scaler  
Scaler



102-242  
Hirschfeld  
file  
3/7 scaler  
anatomical  
hollow  
handle  
6 mm



102-243  
Hirschfeld  
file  
5/11 scaler  
anatomical  
hollow  
handle  
6 mm




102-208  
Colombia  
4R/4L  
Universal  
anatomical  
hollow  
handle  
6 mm




102-209  
Colombia  
13/14  
Universal  
anatomical  
hollow  
handle  
6 mm



102-256  
Black  
26 distel  
scaler  
anatomical  
hollow  
handle  
6 mm



102-1100  
Fischer  
Serration  
anatomical  
hollow  
handle  
6 mm



102-1101  
Pich  
Instruments  
anatomical  
hollow  
handle  
6 mm



102-1102  
Cleoidé  
Discoide  
anatomical  
hollow  
handle  
6 mm



102-256-1  
Black  
27 Mesial  
scaler  
anatomical  
hollow  
handle  
6 mm

102-192-3  
Ostotrimer  
Goldman Fox  
anatomical  
hollow  
handle  
6 mm



*Conservative-Konservierung-Conservador-Soins Conservateur*

---



# AR INSTRUMED

Heidemann Spatulas  
Heidemann-Spatel  
Heidemann-Espátulas  
Spatules de heidemann

Instruments of composite Zirconium Nitride (ZrN) coating tips with anatomical hollow handle 10 mm.

The point covered with (ZrN) ( yellow part) allows to reduce the attachment and the subsequent dragging of the composite in the modeling phase.



102-435

Heidemann Spatula  
1.5mm  
17.5cm



102-436

Heidemann Spatula  
2.0mm  
17.5cm



102-437

Heidemann Spatula  
2.5mm  
17.5cm



102-432

Heidemann Spatula  
2.0mm



102-433

Heidemann Spatula  
2.5mm



102-434

Heidemann Spatula  
3.0mm


Anatomical hollow handle 10 mm

Anatomical hollow handle 6 mm



**AR INSTRUMED**

Burnishers / Filling instruments  
Kugelstopfer / Spatel  
Bruñidores / Atacadores  
Brunissoirs





**AR INSTRUMED**

Burnishers / Filling instruments  
Kugelstopfer / Spatel  
Bruñidores / Atacadores  
Brunissoirs

102-301 Filling Instrument 0.8/1.0 mm Anatomical hollow handle 6 mm	102-302 Filling Instrument 1.2/1.6 mm Anatomical hollow handle 6 mm	102-303 Filling Instrument 2.0/2.5 mm Anatomical hollow handle 6 mm	102-377 Filling Instrument 1.2/1.6 mm Anatomical hollow handle 6 mm	102-394 Hopson 2.0/2.5 mm Anatomical hollow handle 6 mm	102-358 Filling Instrument 1.6/2.5 mm Anatomical hollow handle 6 mm	102-359 Filling Instrument 2.5/3.0 mm Anatomical hollow handle 6 mm	102-399 Mortenson 1.0/1.2 mm Anatomical hollow handle 6 mm
--	--	--	--	--	---	---	---



**AR INSTRUMED**

Curettes Lucas / Hemingway  
Küretten Lucas / Hemingway  
Curetas Lucas / Hemingway  
Curettes de Lucas / Hemingway



103-466

Lucas  
Curette  
2,5mm  
Anatomical  
hollow  
handle 10 mm

102-818

Lucas  
Curette  
2,0mm  
Anatomical  
hollow  
handle 6 mm

102-819  
102-820

Lucas  
Curette  
2,5mm  
Anatomical  
hollow  
handle 6 mm

102-820

Lucas  
Curette  
3,0mm  
Anatomical  
hollow  
handle 6 mm



103-464

Hemingway  
Curette  
2,4mm  
Anatomical  
hollow  
handle 10 mm

102-815

Hemingway  
Curette  
2,5mm  
Anatomical  
hollow  
handle 6 mm

102-816

Hemingway  
Curette  
3,0mm  
Anatomical  
hollow  
handle 6 mm

102-817

Hemingway  
Curette  
3,5mm  
Anatomical  
hollow  
handle 6 mm


Curettes Lucas

Curettes Hemingway



**AR INSTRUMED**

Excavators  
Exkavatoren  
Excavadoras  
Excavateurs




Anatomical hollow handle 10 mm



**AR INSTRUMED**

Excavators  
Exkavatoren  
Excavadoras  
Excavateurs





**AR INSTRUMED**

Cement Spatula  
Zement Spatel  
Espatulas de cemento  
Spatules À ciment



102-484H  
Cement  
Spatula  
6.0mm  
Anatomical  
hollow  
handle 10 mm



102-485H  
Cement  
Spatula  
5.0mm  
Anatomical  
hollow  
handle 10 mm



102-484  
Cement  
Spatula  
6.0mm  
Anatomical  
hollow  
handle 6 mm



102-485  
Cement  
Spatula  
5.0mm  
Anatomical  
hollow  
handle 6 mm



102-486  
Cement  
Spatula  
7.0mm  
Anatomical  
hollow  
handle 6 mm



102-426H  
Spatula  
2.5mm  
Anatomical  
hollow  
handle 10 mm



102-427H  
Spatula  
3.2mm  
Anatomical  
hollow  
handle 10 mm



**AR INSTRUMED**

Filling Instruments  
Füllinstrumente  
Instrumentos para composites  
Instruments à composite



102-395H

Filling Instrument  
1.2 / 1.5mm  
Hollow handle 10 mm



102-396H

Filling Instrument  
1.5 / 2.5mm  
Hollow handle 10 mm



102-312

Filling Instrument  
0.8 / 1.0mm  
Anatomical  
hollow  
handle 6 mm



102-313

Filling Instrument  
1.2 / 1.6mm  
Anatomical  
hollow  
handle 6 mm



102-314

Filling Instrument  
2.0 / 2.5mm  
Anatomical  
hollow  
handle 6 mm



102-409

Filling Instrument  
1.5 / 2.5mm  
Anatomical  
hollow  
handle 6 mm





**AR INSTRUMED**

P.K. Thomas  
P.K. Thomas  
P.K. Thomas  
P.K. Thomas



102-533  
P.K.Thomas  
Anatomical  
hollow handle  
6 mm



102-534  
P.K.Thomas  
Anatomical  
hollow handle  
6 mm



102-535  
P.K.Thomas  
Anatomical  
hollow handle  
6 mm



102-536  
P.K.Thomas  
Anatomical  
hollow handle  
6 mm



102-537  
P.K.Thomas  
Anatomical  
hollow handle  
6 mm




## *Surgery-Chirurgie-Cirugía-Chirurgie*



**AR INSTRUMED**

Scalpel handle  
Skalpell griff  
Mango de bisturí  
Manche de Bistouris



103-004

Scalpel handle  
Long  
170mm

103-001

Scalpel handle  
No.3  
130mm

103-005

Scalpel handle  
No.4  
130mm

106-327

Scalpel handle  
No.3  
with marking 130mm

106-334

Round Scalpel  
handle Adjustable  
170mm

106-325

Round Scalpel  
handle  
Straight

106-326

Round Scalpel  
handle  
Curved

106-336

Ace mini  
Blade handle



# AR INSTRUMED


Premium Line Surgical Scissors

Scissors  
Scheren  
Tijeras  
Ciseaux

High ergonomic design for easy use, thanks to the handle spread apart and bend designed structure.

High hardening production for the serrated and the blade with wire tungsten carbide

The practicality combined with the strength of the instrument allows you to work on particularly hard tissues.





# AR INSTRUMED

## Premium Line Surgical Scissors

Scissors  
Scheren  
Tijeras  
Ciseaux

High ergonomic design for easy use, thanks to the handle spread apart and bend designed structure.

High hardening production for the serrated and the blade with wire tungsten carbide

The practicality combined with the strength of the instrument allows you to work on particularly hard tissues.






**AR INSTRUMED**

*Tungsten Carbide Instruments*

Scissors  
Scheren  
Tijeras  
Ciseaux





**AR INSTRUMED**

*Tungsten Carbide Instruments*

Scissors  
Scheren  
Tijeras  
Ciseaux



103-118-1

Metzenbaum  
Straight in  
Tungsten Carbide  
14.5cm



103-116-1

Metzenbaum  
Curved in  
Tungsten Carbide  
14.5cm



103-113-1

Goldman-Fox  
Straight in  
Tungsten Carbide  
13cm



103-114-1

Goldman-Fox  
Curved in  
Tungsten Carbide  
13cm



**AR INSTRUMED**

*Tungsten Carbide Instruments*

Scissors  
Scheren  
Tijeras  
Ciseaux





**AR INSTRUMED**

Scissors  
Scheren  
Tijeras  
Ciseaux

02-219

Iris  
Straight  
9cm

103-049

Iris  
Straight  
11.5cm

75-095

Iris  
Straight  
13cm

02-221

Iris  
Curved  
9cm

103-050

Iris  
Curved  
11.5cm

75-096


Iris  
Curved  
13cm

103-063

La Grange  
Double Curved  
11.5cm

103-090

La Grange  
Double Curved  
14cm





**AR INSTRUMED**

Scissors  
Scheren  
Tijeras  
Ciseaux



103-110  
Goldman-Fox  
Straight  
13cm



103-112  
Goldman-Fox  
Curved  
13cm



673-2-1  
Goldman-Fox  
Double Curved  
13cm



**AR INSTRUMED**

Scissors  
Scheren  
Tijeras  
Ciseaux




103-065

Dean  
Double Curved  
17 cm

SC-384

Dean  
Double Curved  
14 cm



103-062

Kelly  
Straight  
16 cm



103-107

Kelly  
Curved  
16 cm

103-068

Locklin  
angular  
16 cm

103-056

Martin  
Curved  
11.5 cm



**AR INSTRUMED**

Scissors  
Scheren  
Tijeras  
Ciseaux





**AR INSTRUMED**

Scissors  
Scheren  
Tijeras  
Ciseaux





**AR INSTRUMED**

Scissors  
Scheren  
Tijeras  
Ciseaux



12-002

Lister  
Bandage  
14.5cm

02-009

Dressing  
Scissors  
Straight  
SH/BL  
14cm

02-025

Dressing  
Scissors  
Curved  
SH/BL  
14cm

103-051

Iris  
Scissors  
Angle  
Curved  
11cm

103-101

Universal  
Scissors  
Curved  
16cm



**AR INSTRUMED**

Scissors  
Scheren  
Tijeras  
Ciseaux



103-076-1

Castro Scissors  
Curved  
Tungsten Carbide  
15cm

103-077-1

Castro Scissors  
Curved  
Tungsten Carbide  
18cm

103-074

Westcott  
Scissors  
Curved  
11cm

678-1

Noyes  
Scissors  
Straight  
11cm

678-3

Noyes  
Scissors  
45  
11cm

103-077

Castro  
Scissors  
18cm



**AR INSTRUMED**

*Tungsten Carbide Instruments*

Needle Holders  
Nadelhalters  
Porta-Agujas  
Portes-aiguilles



103-143-1  
Needle holder  
castroviejo  
Barraquer  
straight  
tungsten carbide  
wire 5-0, 6-0,  
7-0  
15cm



25-5007-1  
Needle holder  
castroviejo  
straight  
tungsten carbide  
wire 5-0, 6-0,  
7-0  
14.5cm



25-5006-1  
Needle holder  
castroviejo  
curved  
tungsten carbide  
wire 5-0, 6-0,  
7-0  
14.5cm



25-5009-1  
Needle holder  
castroviejo  
Straight  
tungsten carbide  
wire 5-0, 6-0,  
7-0  
18cm



25-5008-1  
Needle holder  
castroviejo  
curved  
tungsten carbide  
wire 5-0, 6-0,  
7-0  
18cm



**AR INSTRUMED**

*Tungsten Carbide Instruments*

Needle Holders  
Nadelhalter  
Porta-Agujas  
Portes-aiguilles




103-145

Needle holder  
castroviejo  
straight  
tungsten carbide  
wire 5-0, 6-0,  
7-0  
14cm



103-149

Needle holder  
castroviejo  
curved  
tungsten carbide  
wire 5-0, 6-0,  
7-0  
14cm



103-150

Needle holder  
castroviejo  
straight  
tungsten carbide  
wire 5-0, 6-0,  
7-0  
18cm



103-144

Needle holder  
castroviejo  
curved  
tungsten carbide  
wire 5-0, 6-0,  
7-0  
18cm



**AR INSTRUMED**

*Tungsten Carbide Instruments*

Needle Holders  
Nadelhalters  
Porta-Agujas  
Portes-aiguilles





**AR INSTRUMED**

*Tungsten Carbide Instruments*

Needle Holders  
Nadelhalters  
Porta-Agujas  
Portes-aiguilles





**AR INSTRUMED**

*Tungsten Carbide Instruments*

Needle Holders  
Nadelhalters  
Porta-Agujas  
Portes-aiguilles





**AR INSTRUMED**

Needle Holders  
Nadelhalters  
Porta-Agujas  
Portes-aiguilles



11-049

Webster  
Needle Holder  
Tungsten Carbide  
12.5cm

11-116

Gillies  
Needle Holder  
Tungsten Carbide  
16cm

11-070

Ryder  
Needle Holder  
Tungsten Carbide  
13cm

11-070-1

Ryder  
Needle Holder  
Tungsten Carbide  
14cm

11-070-2

Ryder  
Needle Holder  
Tungsten Carbide  
16cm

11-070-3

Ryder  
Needle Holder  
Tungsten Carbide  
18cm

CGF-01

Crown Gripper  
Tungsten Carbide  
16cm



**AR INSTRUMED**

Forceps  
Klemmen  
Pinzas  
Pince



103-173  
Towel Clamp  
Forceps  
11cm



103-174  
Towel Clamp  
Forceps  
13cm




103-170  
Forcep hemostatic  
Straight  
10cm



103-164  
Forcep hemostatic  
Straight  
12cm



103-169  
Forcep hemostatic  
Curved  
10cm



103-165  
Forcep hemostatic  
Curved  
12cm



103-162  
Forcep hemostatic  
Curved  
14cm



102-117  
Peet  
Forceps  
Curved  
11.5cm

05-096  
Mosquito Forceps  
Curved  
1x2  
12cm

05-095  
Mosquito Forceps  
Straight  
1x2  
12cm



**AR INSTRUMED**

Forceps & Needle Holders  
Klemmen & Nadelhalters  
Pinzas & Porta-Agujas  
Pinces & Portes-aiguilles





# AR INSTRUMED

Tweezers  
Pinzetten  
Pinzas  
Pinces



102-104

Surgical tweezers  
ultra fine tips  
straight  
12cm



102-102

Surgical tweezers  
ultra fine tips  
straight  
12cm



102-056  
Tweezer  
for  
tissue  
1:2 Teeth  
15cm



102-086  
De bakey  
tweezer  
atraumatic  
curve for the  
retention of  
mucous tissues  
&  
avoid breakage  
or perforation  
15cm



102-089  
De bakey  
tweezer  
atraumatic  
straight for the  
retention of  
mucous tissues  
&  
avoid breakage  
or perforation  
15cm



102-108  
Surgical tweezers  
ultra fine tips  
curved  
12cm



75-144  
Surgical  
tweezers  
ultra fine tips



**AR INSTRUMED**

Tweezers  
Pinzetten  
Pinzas  
Pinces



102-051-1  
Semken-Taylor  
Straight  
Serrated  
12cm



102-052-1  
Semken-Taylor  
Curved  
Serrated  
12cm



102-053-1  
Semken-Taylor  
Straight  
1:2 Teeth  
12cm



102-054-1  
Semken-Taylor  
Curved  
1:2 Teeth  
12cm



104-478  
Articulating Paper  
Forceps  
Curved  
12cm



102-112  
Articulating Paper  
Forceps  
Straight  
12cm



**AR INSTRUMED**

Tweezers  
Pinzetten  
Pinzas  
Pinces



102-042	102-039	102-041	102-032	102-038	102-040	106-011	04-002	04-004	04-005	102-053	102-056
Adson Forceps Tungsten Carbide 12cm	Adson Micro 12cm	Adson Micro 1x2 12cm	Adson Brown 12cm	Adson Forceps 12cm	Adson Forceps 12cm	Adson Curved 12cm	Dressing Forceps 12cm	Dressing Forceps 14cm	Dressing Forceps 16cm	Tissue Forceps 14cm	Tissue Forceps 15cm



**AR INSTRUMED**

Tweezers  
Pinzetten  
Pinzas  
Pinces



106-025	106-024	106-023	106-026	106-029	106-007	102-110	102-111	103-176	102-068	102-069
Gerald Forceps Straight 18cm	Gerald Forceps Curved 18cm	Gerald Forceps Straight 1x2 18cm	Gerald Forceps Curved 1x2 18cm	Tissue Punch	Suture Plier	Tweezer Wooden Handle Straight 16cm	Tweezer Wooden Handle Curved 16cm	Jone Forceps 9cm	Krane Kaplain 15cm	Krane Kaplain 15cm



**AR INSTRUMED**

Retractors  
Retraktoren  
Separadores  
Ecarteurs



103-197

Senn-Miller  
double ended  
retractor 16cm

103-037

Farabuef  
120x12cm  
(set of 2 pcs)

103-220

Micceldorpf  
retractor  
17x14mm  
21cm

103-221

Micceldorpf  
retractor  
20x22mm  
23cm

103-228

Lengenback  
retractor  
40x11mm  
21cm

103-225

Lengenback  
retractor  
30x11mm  
21cm

103-226

Lengenback  
retractor  
30x14mm  
21cm

103-224

Lengenback  
retractor  
30x16mm  
21cm



**AR INSTRUMED**

Retractors  
Retraktoren  
Separadores  
Ecarteurs

103-183

Sternberg  
retractor  
16cm



103-182

Sternberg  
retractor  
14cm

104-856

Retractor  
Universal  
14cm

103-185

Orringer  
retractor  
Medium  
11x8cm

103-024

Stainless steel  
double-ended  
cheek retractor

104-591

Bishop retractor

103-189

Brunings retractor  
19cm

106-248

Minnesota  
retractor  
large



**AR INSTRUMED**

Retractors/Matrix Retainers  
Retraktoren/Matrizenspanner  
Separadores/Matrix Retainers  
bouche-Gag /Matrice type tofflemire



106-238

Molt  
mouth gag  
11cm

106-236

Molt  
mouth gag  
14cm



103-368  
Tofflemire  
universal  
Fig.1

103-369  
Tofflemire  
Junior  
Fig.2

103-370  
Tofflemire  
Senior  
Fig.3

103-364  
Ivory 8  
Fig.6  
8cm

103-364-1  
Nystrom  
Right

103-364-2  
Nystrom  
Left

103-383  
Ivory 1



**AR INSTRUMED**

Raspatories  
Raspatores  
Separadores  
Pistolets/Serringues à Amalgame



103-438

Williger  
13cm



103-437

Williger  
16cm



103-435

Williger  
17cm

103-436

Williger  
17cm

103-451


Williger  
bone  
curette  
1/2

103-450

Williger  
bone  
curette  
0/1

103-449

Williger  
bone  
curette  
00/0



103-352

Amalgam  
carrier  
14cm 90°

103-351

Amalgam  
carrier  
14cm 45°

103-350

Amalgam  
carrier  
14cm 15°



2.0mm

103-339

Amalgam  
carrier  
20cm



2.0mm

103-340

Amalgam  
carrier  
20cm



**AR INSTRUMED**


Rubber dam Forceps/Punchs  
Kofferdam klammerzange/Kofferdamlochzange  
Porta Clamps Pliers/Perforador de Diques  
Pince Porte Clamp/Pince à Digue





**AR INSTRUMED**

Rubber dam Clamps  
Kofferdam  
Porta Clamps  
Set de Pince Porte Clamp/Pinces à Digue



103-267-1  
Rubber dam Clamps  
Set of 13 Pcs  
1 Tray



103-299 103-301 103-303 103-304

Clamp No 4 Clamp No 7 Clamp No 8 Clamp No 8a

103-308 103-309 103-310 103-311

Clamp No 12a Clamp No 13a Clamp No 14 Clamp No 14a

103-304 103-308 103-308 103-284-1

Clamp No 27 Clamp No 2 Clamp No 00 Clamp No 26N



103-486  
X-Ray Holder  
Clips Single



102-1000  
Finger Protector  
Small

102-1001  
Finger Protector  
Medium

102-1002  
Finger Protector  
Large

103-488  
X-Ray Holder  
Set of 10  
Clips

103-601  
Napkin Holder  
Stainless Steel



103-268  
Rubber dam frame  
for adults  
9x9cm

103-267  
Rubber dam frame  
for adults  
9x10cm




103-362  
Amagam Well  
Stainless Steel



**AR INSTRUMED**

Retractors  
Retraktoren  
Separadores  
bouche-Gag





*Micro Surgery-MikroChirurgie-Micro Cirugía-Micro Chirurgie*



**AR INSTRUMED**

Micro Needle Holders/Scissors  
Mikro Nadelhalter/Scheren  
Micro Porta-Agujas/Tijeras  
Micro Portes-aiguilles / Ciseaux



K25-5007  
Micro Castroviejo  
Needle Holder  
Tungsten Carbide  
str  
140mm



K25-5006  
Micro Castroviejo  
Needle Holder  
Tungsten Carbide  
cvd  
140mm



K25-5008  
Micro Castroviejo  
Needle Holder  
Tungsten Carbide  
cvd  
180mm



109-005  
Micro Castroviejo  
Scissors  
Tungsten Carbide  
cvd  
180mm



**AR INSTRUMED**

Micro Scissors  
Mikro Scheren  
Micro Tijeras  
Micro Ciseaux





**AR INSTRUMED**

Micro Scissors  
Mikro Scheren  
Micro Tijeras  
Micro Ciseaux





#### Precleaning and Cleaning

Removing bioburden and any adherent impression material is the fundamental first step in achieving an effective disinfection and sterilization.

#### Disinfecting

Resistant to disinfecting solutions with high quality AISI 304 high stainless steel.

#### Sterilizing

Autoclavable at 134°C (273,2°F) with no alteration.

## *Impression Trays - Abdrucklöffel - Cubetas de Impresión-Portes empreintes*



**AR INSTRUMED**

Perforated/Solid  
Perforierte/Solide  
Perforada/Solido  
Perfore/Plein



**Perforated Lower Jaw**

- |        |        |         |
|--------|--------|---------|
| Fig. 1 | Small  | 107-351 |
| Fig. 2 | Medium | 107-352 |
| Fig. 3 | Large  | 107-353 |



**Perforated Upper Jaw**

- |        |        |         |
|--------|--------|---------|
| Fig. 1 | Small  | 107-361 |
| Fig. 2 | Medium | 107-362 |
| Fig. 3 | Large  | 107-363 |



**Solid Lower Jaw**

- |        |        |         |
|--------|--------|---------|
| Fig. 1 | Small  | 107-371 |
| Fig. 2 | Medium | 107-372 |
| Fig. 3 | Large  | 107-373 |



**Solid Upper Jaw**

- |        |        |         |
|--------|--------|---------|
| Fig. 1 | Small  | 107-381 |
| Fig. 2 | Medium | 107-382 |
| Fig. 3 | Large  | 107-383 |



**AR INSTRUMED**

Perforated/Solid  
Perforierte/Solide  
Perforada/Solido  
Perfore/Plein



**Perforated Lower Jaw**

- |        |        |         |
|--------|--------|---------|
| Fig. 1 | Small  | 107-391 |
| Fig. 2 | Medium | 107-392 |
| Fig. 3 | Large  | 107-393 |



**Perforated Upper Jaw**

- |        |        |         |
|--------|--------|---------|
| Fig. 1 | Small  | 107-401 |
| Fig. 2 | Medium | 107-402 |
| Fig. 3 | Large  | 107-403 |

**Solid Lower Jaw**

- |        |        |         |
|--------|--------|---------|
| Fig. 1 | Small  | 107-411 |
| Fig. 2 | Medium | 107-412 |
| Fig. 3 | Large  | 107-413 |

**Solid Upper Jaw**

- |        |        |         |
|--------|--------|---------|
| Fig. 1 | Small  | 107-421 |
| Fig. 2 | Medium | 107-422 |
| Fig. 3 | Large  | 107-423 |



**AR INSTRUMED**

Perforated/Solid  
Perforierte/Solide  
Perforada/Solido  
Perfore/Plein



**Solid Upper Jaw**

- Fig. 1 Small 107-040  
Fig. 2 Medium 107-039  
Fig. 3 Large 107-038  
Fig. 4 Extra Large 107-037

**Solid Lower Jaw**

- Fig. 1 Small 107-045  
Fig. 2 Medium 107-044  
Fig. 3 Large 107-043  
Fig. 4 Extra Large 107-042

**Rimlock**



**Perforated Upper Jaw**

- Fig. 1 Small 107-004  
Fig. 2 Medium 107-003  
Fig. 3 Large 107-002  
Fig. 4 Extra Large 107-001

**Perforated Lower Jaw**

- Fig. 1 Small 107-010  
Fig. 2 Medium 107-009  
Fig. 3 Large 107-008  
Fig. 4 Extra Large 107-007



**AR INSTRUMED**

Perforated/Solid  
Perforierte/Solide  
Perforada/Solido  
Perfore/Plein



**Solid Lower Jaw**

- Fig. 1 Small 107-113  
Fig. 2 Medium 107-112  
Fig. 3 Large 107-111  
Fig. 4 Extra Large 107-110

**Rimlock**



**Perforated Upper Jaw**

- Fig. 1 Small 107-095  
Fig. 2 Medium 107-094  
Fig. 3 Large 107-093  
Fig. 4 Extra Large 107-092

**Perforated**

L2D	107-016	107-013
L4D	107-017	107-014
L6D	107-018	107-015



## Implantology-Implantologie-Implantología-Implantologie



# AR INSTRUMED

Bone scraper / Replacement blade / Application Syringe  
Knochenschaber / Ersatzklinge fur / Applikations Spritze  
Raspadores de hueso / Cuchilla de respuesto / Jeringa de Aplicación  
Os-Grattoir/ Seringue À os



H2020  
Compactor spoon for  
intake and bone  
Ø 6.1mm  
probe  
with marks



106-339  
Bone Scraper  
Curved  
15.5cm  
  
106-339-2  
Replacement blade  
of curved  
bone scraper

*Bone Compactor*

*Bone Scraper*



*Application Syringe*



**AR INSTRUMED**

Tissue Punches  
Hautstanzen  
Pince à couper le tissu  
Poinçons de tissus



106-290-1

Tissue Punch  
Straight  
3mm

106-290-2

Tissue Punch  
cvd  
3mm

106-290-3

Tissue Punch  
Straight  
4mm

106-290-4

Tissue Punch  
cvd  
4mm

106-290-5

Tissue Punch  
Straight  
5mm

106-290-6

Tissue Punch  
cvd  
5mm

106-290-7

Tissue Punch  
Straight  
6mm

106-290-8

Tissue Punch  
cvd  
6mm



# AR INSTRUMED

Implant Instrumets  
Implantologie Instrumente  
Instrumentos de implantologia  
Implantologie Instruments



106-272

Bone Morsellezr  
Titanium

Bone mill with replaceable Titanium cutting surface  
for preparation of autogenous bone.

Knochenmühle mit auswechselbarer Titan Zähnung zur  
gewinnung von feinporen autogem kochen.


Molino de hueso con diente de Titano reemplazables  
para la reparicon de hueso autologo.



75-219

Bone Morsellezr  
in Titanium  
Ø 40mm

Bone cracker , for longer bone pieces  
knochenschlagdose , fur längere knochenstücke  
Triturador de hueso , para piezas de hueso mas grande



106-275


Mixing bowl  
for bone material  
Mischdose fur  
knochenmaterial  
Caja de mezcla para  
material de hueso





**AR INSTRUMED**

Curettes Gracey  
Gracey Küretten  
Curetas Gracey  
Curettes de gracey



Curettes Gracey with anodization of Titanium (Blue) coating tips with anatomical hollow handle 10 mm



**AR INSTRUMED**

Periosteal Elevator  
Raspatorium  
Periostotomos  
Périotomes



103-424

Buser  
17,5cm

106-049

Buser  
18cm

405

Double Pichard  
18,5cm

106-051H

Prichard  
19cm




Anatomical hollow handle 10 mm



**AR INSTRUMED**

Periosteal Elevator  
Raspatorium  
Periostotomos  
Périotomes



Anatomical hollow handle 10 mm



**AR INSTRUMED**

Osteotomy Set  
Osteotomie Set  
Set de osteotomía  
Set de Ostéotomes



106-465

106-466

106-467

106-468

106-469

106-470

106-471

106-472

106-473

106-474

Osteotome  
Straight-Concave  
Marks:  
1,0-1,2-1,4-1,6-1,8  
Ø 4,9 -180mm

Osteotome  
Straight-Concave  
Marks:  
1,0-1,2-1,4-1,6-1,8  
Ø 4,1 -180mm

Osteotome  
Straight-Concave  
Marks:  
1,0-1,2-1,4-1,6-1,8  
Ø 3,2 -180mm

Osteotome  
Straight-Concave  
Marks:  
1,0-1,2-1,4-1,6-1,8  
Ø 2,7 -180mm

Osteotome  
Straight-Concave  
Marks:  
1,0-1,2-1,4-1,6-1,8  
Ø 2,0 -180mm

Osteotome  
Curved-Concave  
Marks:  
1,0-1,2-1,4-1,6-1,8  
Ø 4,9 -180mm

Osteotome  
Curved-Concave  
Marks:  
1,0-1,2-1,4-1,6-1,8  
Ø 4,1 -180mm

Osteotome  
Curved-Concave  
Marks:  
1,0-1,2-1,4-1,6-1,8  
Ø 3,2 -180mm

Osteotome  
Curved-Concave  
Marks:  
1,0-1,2-1,4-1,6-1,8  
Ø 2,7 -180mm

Osteotome  
Curved-Concave  
Marks:  
1,0-1,2-1,4-1,6-1,8  
Ø 2,0 -180mm



**AR INSTRUMED**

Osteotomy Set  
Osteotomie Set  
Set de osteotomía  
Set de Ostéotomes



106-475	106-476	106-477	106-478	106-479	106-480	106-481	106-482	106-483	106-484
Osteotome Straight-Convex Marks: 1,0-1,2-1,4-1,6-1,8 Ø 4,9 -180mm	Osteotome Straight-Convex Marks: 1,0-1,2-1,4-1,6-1,8 Ø 4,1 -180mm	Osteotome Straight-Convex Marks: 1,0-1,2-1,4-1,6-1,8 Ø 3,2 -180mm	Osteotome Straight-Convex Marks: 1,0-1,2-1,4-1,6-1,8 Ø 2,7 -180mm	Osteotome Straight-Convex Marks: 1,0-1,2-1,4-1,6-1,8 Ø 2,0 -180mm	Osteotome Curved-Convex Marks: 1,0-1,2-1,4-1,6-1,8 Ø 4,9 -180mm	Osteotome Curved-Convex Marks: 1,0-1,2-1,4-1,6-1,8 Ø 4,1 -180mm	Osteotome Curved-Convex Marks: 1,0-1,2-1,4-1,6-1,8 Ø 3,2 -180mm	Osteotome Curved-Convex Marks: 1,0-1,2-1,4-1,6-1,8 Ø 2,7 -180mm	Osteotome Curved-Convex Marks: 1,0-1,2-1,4-1,6-1,8 Ø 2,0 -180mm



# AR INSTRUMED

Sinus Lift Instruments  
Sinus Lift Instrumente  
Instrumentos Sinus Lift  
Kit pour sinus lift



106-059-1H

Sinus Lift  
Instrument  
Fig. 1

106-059-2H

Sinus Lift  
Instrument  
Fig. 2

106-059-3H

Sinus Lift  
Instrument  
Fig. 3


106-059-4H

Sinus Lift  
Instrument  
Fig. 4



Sinus floor elevators, tightly curved, for raising the deep area of the membrane from the floor of the sinus cavity.

Sinus lift instruments anodization of Titanium (Blue) coating tips with anatomical hollow handle 10 mm



106-117

Molt 4  
With Suture  
Hole  
10mm Hollow Handle

103-429

Freer  
10mm Hollow  
Handle

104-495

Bone Graft  
3.3/4mm  
10mm Hollow  
Handle

104-496

Graft/Carrier  
2-4-6-8-10  
10mm Hollow  
Handle



**AR INSTRUMED**


Bone Scraper / Chisel  
Knochenschaber / Meißel  
Raspador de hueso / Cinceles  
Burin dentaire

106-340 Chisel 1mm / 4mm	106-362 Bone Chisel Straight Centre 5mm	106-363 Bone Chisel Left Curved 5mm	106-364 Bone Chisel Right Curved 5mm	106-339-1 Bone Scraper <p>The Bone Scraper has extended working area for bone curettage so that a large volume of bone substance can be obtain. The shortened curved shaft ideally meets the surgeons needs in handling the jaw under different anatomical conditions.</p>	35-1250 Tweezer for implant 18cm	106-131 Seldin 18cm
--------------------------------	---	---	--	--	--	---------------------------




**AR INSTRUMED**

Suction pipe/Compactor  
Chirurgischer absauger/Stopfer  
Tubo de succión/compacteurs  
Compacteur À os



100-322  
Suction tube  
(Titanium tube)  
 $\varnothing$  2mm  
15cm



100-315  
Suction pipe  
Small  
 $\varnothing$  3mm  
18cm



100-316  
Suction pipe  
Medium  
 $\varnothing$  4mm  
18cm



100-317  
Suction pipe  
Large  
 $\varnothing$  5mm  
18cm



106-285  
Compactor  
Serrated  
 $\varnothing$  3mm  
 $\varnothing$  4mm



106-284  
Compactor  
Serrated  
 $\varnothing$  6mm  
 $\varnothing$  8mm



106-283  
Compactor  
Bone Spoon  
 $\varnothing$  3mm  
 $\varnothing$  6mm




106-214H  
Bone  
File

Anatomical hollow handle 10 mm



**AR INSTRUMED**

Hammer  
Mallet  
Martillo  
Maillets




**106-310**  
Tissue nipper  
13cm



**AR INSTRUMED**

Implant Instruments  
Implantologie Instrumente  
Instrumentos de implantología  
Implantologie Instruments





## *Orthodontics-Kieferorthopädie-Ortodoncia-Orthodontie*

---



## AR INSTRUMED

This instrument cuts wires close to the buccal tube and safely holds the distal end. The Tungsten Carbide inserts ensure precision cutting, long wear and unsurpassed cutting performance. Cutting Capacity: archwires up to .021" x .025" (.55 x .64mm)



This instrument cuts wires close to the buccal tube and safely holds the distal end. The Tungsten Carbide inserts ensure precision cutting, long wear and unsurpassed cutting performance. Cutting Capacity: archwires up to .012" x .020" (.30 x .50mm)





# AR INSTRUMED

Orthodontics Pliers  
Kieferorthopädie Zangen  
Alicates de Ortodoncia  
Pince coupante d'orthodontie

Standard cutter with 7 degree angled tips to allow easy access and better visibility in difficult to reach areas. Inserted tips ensure precision cutting of ligature wire, pins and elastics.  
Maximum Cutting Capacity: soft wire up to .016" (.41mm)



105-159-1

Pin & Ligature  
Soft Wire Cutter  
Tungsten-carbide

Easily cuts any size archwire or pin. The inserted tips ensure precision cutting and perfect tip alignment.  
Minimum Cutting Capacity: .012" (.30mm)  
Maximum Cutting Capacity: archwires up to .021" x .025" (.55mm x .64mm)



105-159

Hard wire cutter  
Plier  
Tungsten-carbide



## AR INSTRUMED

The most popular wire bending plier. The inserted tips provide excellent wear and bend archwires from .012" to .020" (.30mm – .51mm). Beaks are grooved for round wire at 1mm and 2mm from tip.  
Cutter at the base of beaks cuts archwires from .014" to .020" (.36mm – .51mm).



105-090

Light Wire bending  
Plier with Cutter  
Tungsten-carbide



105-162

Heavy Duty Hard  
Wire Cutter  
Tungsten-carbide

Orthodontics Pliers  
Kieferorthopädie Zangen  
Alicates de Ortodoncia  
Pince coupante d'orthodontie

This instrument cuts up to 1.5mm wire with ease due to its triple action design.  
Inserted tips ensure precision cutting and perfect tip alignment.

Minimum Cutting Capacity: .014" (.36mm)

Maximum Cutting Capacity: wire up to .060" (1.5mm)



**AR INSTRUMED**

Orthodontics Pliers  
Kieferorthopädie Zangen  
Alicates de Ortodoncia  
Pince coupante d'orthodontie

Precision serrated tips are designed for heavy wires and to prevent the wire from slipping. For safe and easy insertion or removal of all arches. Holds wire up to .030" (.76mm).





# AR INSTRUMED

Posterior Band Remover with a longer lower band removing arm.

Replacement pads are available in:  
3/16" (5mm) and 1/8" (4mm) diameters.



Orthodontics Pliers  
Kieferorthopädie Zangen  
Alicates de Ortodoncia  
Pince coupante d'orthodontie

The serrated tips of this instrument hold the wire securely while bending the wire 90 degree distal to the buccal tube, avoiding patient discomfort.

For use with Nickel Titanium wires, .012" – .018" (.30 x 46mm).





# AR INSTRUMED

Orthodontics Pliers  
Kieferorthopädie Zangen  
Alicates de Ortodoncia  
Pince coupante d'orthodontie

Force is applied to the bracket tooth interface. The tips are positioned mesially, distally or at the occlusal gingival interface to remove direct bond brackets safely and comfortably.



Direct Bond Remover (Angled) with 60 degree angle allows easier access in those difficult to reach posterior areas.



# AR INSTRUMED

Precision aligned tips are rounded for consistent contouring and bending without damaging the archwire. One-piece construction of the double tip ensures superior strength. Can be used on wire up to .030" (.76mm).



Orthodontics Pliers  
Kieferorthopädie Zangen  
Alicates de Ortodoncia  
Pince coupante d'orthodontie

Designed for surgical ball hooks and auxiliaries.  
The slotted tip holds auxiliaries secure to provide for positive crimping





## AR INSTRUMED

Precision loop forming with four step tips 3mm, 4mm, 5mm and 6mm.  
Excellent for round and square wires up to .022" (.56mm).  
All working edges have been bevelled to prevent wire scoring.



Orthodontics Pliers  
Kieferorthopädie Zangen  
Alicates de Ortodoncia  
Pince coupante d'orthodontie

Designed for torquing and bending square or rectangular wires. The blades are designed to be parallel at .020" opening.  
Can be used on wires up to .021" x .025" (.55 x .64mm).





# AR INSTRUMED

Precision matched concave and convex tips for smooth contouring of light wire loops. Contours and forms arches without torquing.  
Can be used on wires up to .021" x .025" (.55mm x .64mm).



105-098  
De la rosa  
Plier



105-216  
Young  
Plier

Orthodontics Pliers  
Kieferorthopädie Zangen  
Alicates de Ortodoncia  
Pince coupante d'orthodontie


For labial, lingual, arches and canine loops.  
For hard wire up to .028" (0.7mm) max.



**AR INSTRUMED**

Orthodontics Pliers  
Kieferorthopädie Zangen  
Alicates de Ortodoncia  
Pince coupante d'orthodontie

Wire bending and loop forming plier  
Maximum Bend wire up to: .028" (0.70mm)



Wire bending and loop forming plier  
Maximum Bend wire up to: .028" (0.70mm)





# AR INSTRUMED

The most popular and versatile pliers for working with round wire up to .30" (.76mm) in diameter.

The round beak is precision ground to a perfect cone, .025" (.64mm) at the tip.  
The beaks are parallel at .020" (.51mm) opening.



Longer beaks for making precision loops in light wires.  
Can be used with round and square wire up to .020" (.51mm).



Orthodontics Pliers  
Kieferorthopädie Zangen  
Alicates de Ortodoncia  
Pince coupante d'orthodontie



**AR INSTRUMED**

Orthodontics Pliers  
Kieferorthopädie Zangen  
Alicates de Ortodoncia  
Pince coupante d'orthodontie





# AR INSTRUMED

Micro-miniature light wire cutter, it cuts soft ligature wires, for use in case of limited access. Caution Maximum Cutting Capacity: 0.3mm (.012")



105-159-2  
Soft Wire cutter  
Plier

Universal dolphin style plier for bending wires up to 0.9mm (.035") or cutting up to 0.7mm (.028")



105-156  
Universal Dolphin  
Style Plier



105-203  
Schwarz  
Arrow Clasp  
Forming & Bending  
Pliers

Orthodontics Pliers  
Kieferorthopädie Zangen  
Alicates de Ortodoncia  
Pince coupante d'orthodontie



**AR INSTRUMED**

Orthodontics Pliers  
Kieferorthopädie Zangen  
Alicates de Ortodoncia  
Pince coupante d'orthodontie



105-170  
Johnson  
Pliers



105-092  
How  
Utility  
Plier





**AR INSTRUMED**

Orthodontics Pliers  
Kieferorthopädie Zangen  
Alicates de Ortodoncia  
Pince coupante d'orthodontie



105-210  
Wire Bending  
Loop Forming  
Pliers



105-215  
Schwarz  
Medium  
Forming & Bending  
Plier



**AR INSTRUMED**

Orthodontics Pliers  
Kieferorthopädie Zangen  
Alicates de Ortodoncia  
Pince coupante d'orthodontie





**AR INSTRUMED**

Orthodontics Pliers  
Kieferorthopädie Zangen  
Alicates de Ortodoncia  
Pince coupante d'orthodontie



105-269-1

Micro  
Wire Cutter  
Tungstan Carbide  
Plier  
Long Handle



105-266

Distal End Cutter  
Tungstan Carbide  
Plier  
Long Handle



105-267

Distal End Cutter  
Curved Handle  
Tungstan Carbide  
Plier  
Long Handle



**AR INSTRUMED**

Orthodontics Pliers  
Kieferorthopädie Zangen  
Alicates de Ortodoncia  
Pince coupante d'orthodontie



105-273  
Posterior Band  
Remover  
Plier  
Long Handle



105-271  
Weingart  
Plier  
Long Handle



105-270  
Direct Bond  
Remover Angled  
Plier  
Long Handle



**AR INSTRUMED**

Orthodontics Pliers  
Kieferorthopädie Zangen  
Alicates de Ortodoncia  
Pince coupante d'orthodontie





**AR INSTRUMED**

Ligature Forcep / Applicator for Elastic Ligature  
Ligaturspannzange / Elastik-Modul Adapter  
Pinza para Ligaduras / Aplicador de Ligaduras Elasticas  
Instruments Orthodontie

Pliers and applicators for elastic





**AR INSTRUMED**

Ligature Forcep / Applicator for Elastic Ligature  
Ligaturspannzange / Elastik-Modul Adapter  
Pinza para Ligaduras / Aplicador de Ligaduras Elasicas  
Instruments Orthodontie

### Ligature Instruments



105-060  
Directional  
for Ligature  
17.5cm




105-051  
Instrument  
utility Ligature  
17.5cm



75-202  
Direct bonding  
bracket  
holder  
14cm




105-014  
Bracket  
positioning  
tweezer  
12cm



105-021  
Measurer for bands  
and brackets  
with  
different heights



105-002  
Easy access  
tweezer  
12cm



105-007  
Delux bracket  
holder  
12cm



105-001  
Bracket  
Placing  
Tweeze  
13cm



105-012  
Bracket  
Tweezer  
12cm



**AR INSTRUMED**

Orthodontics Pliers  
Kieferorthopädie Zangen  
Alicates de Ortodoncia  
Pince coupante d'orthodontie



105-022

Bracket  
positioning  
gauge  
.018" slot

105-023


Bracket  
positioning  
gauge  
.022" slot



105-501

Stand for Pliers  
Stainless Steel





*Prosthesis-Prothèse-Prótesis-Prthèse*

---



# AR INSTRUMED

Forceps & Pliers  
Zahnzangen  
Fórceps para extracción  
Daviers & Pinces

Precision pliers for removal of temporary crowns and ceramic bridges, can be used both for the lower jaw that inserts for use with interchangeable part to be removed completely dry





**AR INSTRUMED**

Forceps & Pliers  
Zahnzangen  
Fórceps para extracción  
Daviers & Pinces


These pliers allow for careful, precise spreading of the split crown.  
Upon spreading the split, the crown can be easily removed.



100-331  
Crown Remover  
Forceps  
Left



100-332  
Crown Remover  
Forceps  
Right



103-253  
Crown spreader  
plier



**AR INSTRUMED**

Crown Remover  
Kronenentferner  
Martillo Arranca Coronas  
Arrache Couronne

Plier for recovery fragments of root canal instruments and  
for inserting fiber posts.



100-237  
Pliers for  
fragments  
of root canal  
14cm



103-238  
Bohm  
11.5cm



103-236  
Furrer  
Copper ring  
removing forceps  
16cm  
Straight



**AR INSTRUMED**

Crown Remover  
Kronenentferner  
Martillo Arranca Coronas  
Arrache Couronne



103-247

Crown Remover  
Complete with 3 Tips and Mounting kit



103-248

Crown Remover Automatic  
Complete with 3 Tips and Mounting kit



103-249


Crown Remover  
Complete with 5 Attachments



**AR INSTRUMED**

Crown Remover  
Kronenentferner  
Martillo Arranca Coronas  
Arrache Couronne





## *Calipers-Messinstrumente-Calibres-Instruments de mesure*


---



# AR INSTRUMED

Calipers  
Messzirkel  
Calibre  
Instruments de mesure


Caliper graduated in MM  
Schieblehre, graduiert in MM  
Calibrador, graduado en MM





**AR INSTRUMED**

Calipers  
Messzirkel  
Calibre  
Instruments de mesure



106-216  
Calipers  
0 - 90mm

Mini dental caliper for  
width and depth mapping



106-216-1  
Caliper  
0 - 80mm



106-221  
Bone  
calliper  
gauge  
35mm



104-998  
Bone  
calliper  
gauge

Ridge mapping caliper  
for mapping bone width  
Messzange (Osteometer)  
zur kieferkammbreiten  
messung  
Calibre de implantología  
para medir la anchura  
de hueso



## *Laboratory-Laboratorium-Laboratorio-Laboratoire*



**AR INSTRUMED**

Spatulas / Knives  
Spatcheln / Messer  
Espátulas / Cuchillas  
Spatules et couteaux



102-994

Fahnenstock  
13cm

102-995

Fahnenstock  
17cm

102-998

Lessmann  
13cm

102-999

Lessmann  
17cm

102-996

Knife  
12cm

102-997

Knife  
17cm

102-1005

Gritmann  
16cm

102-1006

Knife  
18cm



**AR INSTRUMED**

Spatulas / Knives  
Spackeln / Messer  
Espátulas / Cuchillas  
Spatules et couteaux



102-1014

Spatula  
flexible  
19.5cm



102-1014-1

Spatula  
flexible  
19.5cm



102-1016

Spatula  
flexible  
21.5cm



102-1008

Spatula  
very flexible  
19cm



102-1011

Spatula  
lichtenstein  
17cm



103-482

Plaster cutting  
scissor  
17cm



# AR INSTRUMED

Knives  
Messer  
Cuchillas  
couteaux



102-520  
Hylin  
le corn  
anatomical  
hollow  
handle  
6 mm

102-521  
Hylin  
le corn  
anatomical  
hollow  
handle  
6 mm

102-522  
Craving  
instrument  
anatomical  
hollow  
handle  
6 mm

102-517  
Zahle  
anatomical  
hollow  
handle  
6 mm

102-518  
Zahle  
anatomical  
hollow  
handle  
6 mm

102-520-1  
Craving  
instrument  
anatomical  
Solid  
handle  
6 mm

102-520-2  
Craving  
instrument  
anatomical  
Solid  
handle  
6 mm

102-520-3  
Craving  
instrument  
anatomical  
Solid  
handle  
6 mm

102-520-4  
Craving  
instrument  
anatomical  
Solid  
handle  
6 mm

102-520-5  
Craving  
instrument  
anatomical  
Solid  
handle  
6 mm

102-520-6  
Craving  
instrument  
anatomical  
Solid  
handle  
6 mm



*Sterilization-Sterilisation-Esterilización-Stérilisation*

---



# AR INSTRUMED

Sterilization trays of stainless steel with locking system  
& silicone autoclavable.

107-147

10 Pcs.  
Sterilization tray  
182x144x34mm

107-145


5 Pcs.  
Sterilization tray  
182x72x34mm

107-146

20 Pcs.  
Sterilization tray  
288x182x34mm

107-141

10 Pcs.  
Sterilization tray  
140x203x32mm






**AR INSTRUMED**

Sterilization Trays  
Kassetten Halten  
Esterilización Bandejas  
Plateaux de stérilisation



107-150-1

5Pcs  
Sterilization tray  
182x72x34cm



107-151-1

10Pcs  
Sterilization tray  
182x144x34cm



**AR INSTRUMED**

Instruments Tray/Box  
Instrumenten Tablet/Box  
Bandejas Cajas Instrumentos  
Instruments Barquettes & Boîtes



107-193

Instrument tray  
20x10x2cm



107-193-1

Instrument tray  
20x15x2cm



107-241

Instrument tray  
25cm



107-339

Endodontic box  
round



107-340

Endodontic box  
Rectangular 10x4x5cm



107-205

Instrument box  
Rectangular 20x10x5cm



**AR INSTRUMED**

Stainless Steel Bowl  
Edelstahl Schüssel  
Acero inoxidable Tazón  
Bol en acier inoxydable

Containers for mixing.





*Dental Implant - Zahnimplantat - Implante dental - Implant dentaire*

---



**AR INSTRUMED**

Dental Implant  
Zahnimplantat  
Implante dental  
Implant dentaire

## Bone Expander Kit



108-001

Bone Expander  
Kit  
Regular

Dental Implant Bone Expander Sinus Lift Kit, a Best Alternative To Osteotomes For The Expansion And Condensing Of The Atrophic Mandible And Maxilla.

**The Kit Includes:**

- 5 x Bone Expansion Screws (2.6mm, 3.0mm, 3.4mm, 3.8mm, 4.3mm)
- 1 x Ratchet Extender
- 1 x 4.0mm Motor Mount Adapter
- 1 x Hand Driver
- 1 x 1.8mm Pilot Drill
- 2 x Parallel Pins
- 1 x Ratchet Wrench, 4.0mm Square Driver
- 1 x Sterilization Cassette / Box



108-002

Bone Expander  
Kit

Dental Implant Bone Expander Sinus Lift Kit, a Best Alternative To Osteotomes For The Expansion And Condensing Of The Atrophic Mandible And Maxilla.

**The Kit Includes:**

- 5 x Bone Expansion Screws (2.6mm, 3.0mm, 3.4mm, 3.8mm, 4.3mm)
- 1 x Ratchet Extender 4.0mm
- 1 x 4.0mm Motor Mount Adapter (for Handpiece)
- 1 x Hand Driver / Wrench
- 1 x 1.6mm Pilot Drill
- 2 x Saw Disks (7.0mm & 10.0mm)
- 1 x Ratchet Wrench, 4.0mm Square Driver (Long)
- 1 x Sterilization Cassette / Box



**AR INSTRUMED**

Dental Implant  
Zahnimplantat  
Implante dental  
Implant dentaire

## Bone Compression Kit



108-003

### Bone Compression Kit

The Bone Compression Kit enhances primary Stability in Type 3 or 4 bone by gently condensing the bone prior to implant placement.

The unique design of the compression screws also facilitates simple and effective closed sinus augmentation procedures.

#### Special Advantages:

- Prevents traumatic Osteotomy
- Increase bone density
- Increases Primary Stability



108-004

### Bone Compression Kit



**AR INSTRUMED**

Dental Implant  
Zahnimplantat  
Implante dental  
Implant dentaire

## Basic Dental Implant Kit



108-005

Basic Dental Implant  
Kit

Basic Dental Implant Universal Kit, consisting of 25pcs Essential Tools & Instruments for an Implant Surgery.

**The Kit Includes:**

- |                           |                             |
|---------------------------|-----------------------------|
| 8pcs Trehpine Drills      | 8pcs Titanium-Coated Drills |
| - Int: 2.0mm, Ext: 3.0mm  | - Dia: 2.0mm, L: 18.5mm     |
| - Int: 3.0mm, Ext: 4.0mm  | - Dia: 2.5mm, L: 18.5mm     |
| - Int: 4.0mm, Ext: 5.0mm  | - Dia: 2.8mm, L: 18.5mm     |
| - Int: 5.0mm, Ext: 6.0mm  | - Dia: 3.2mm, L: 18.5mm     |
| - Int: 6.0mm, Ext: 7.0mm  | - Dia: 3.65mm, L: 18.5mm    |
| - Int: 7.0mm, Ext: 8.0mm  | - Dia: 4.5mm, L: 18.5mm     |
| - Int: 8.0mm, Ext: 9.0mm  | - Dia: 4.8mm, L: 18.5mm     |
| - Int: 9.0mm, Ext: 10.0mm | - Dia: 5.2mm, L: 18.5mm     |
- Drill Extender
  - Driver 2.5mm
  - Compression Screw 2.8mm
  - Compression Screw 4.0mm
  - Ratchet Driver 1.25 mm
  - Ratchet Driver 2.5mm
  - Hand Driver 1.25mm
  - Hand Driver 2.5mm
  - Torque Wrench (10-40Ncm)



108-006

Basic Dental Implant  
Kit

Consisting 50pcs Dental Implant Tools & Instruments with Special Conical Drills including Stoppers, this is a very useful Kit for an Implant Doctor.

**16pcs Conical Drills with Stoppers**

- 1.8mm / 2.2mm - S1
- 1.8mm / 2.2mm - S3
- 1.8mm / 2.2mm - S5
- 1.8mm / 2.2mm - S7
- 2.2mm / 2.7mm - S1
- 2.2mm / 2.7mm - S3
- 2.2mm / 2.7mm - S5
- 2.2mm / 2.7mm - S7

**Bone Expander Screws**

- 1.6mm / 2.6mm
- 1.9mm / 3.0mm
- 2.9mm / 3.4mm
- 2.3mm / 3.8mm
- 2.9mm / 4.3mm
- 2.8mm / 4.5mm - S1
- 2.8mm / 4.5mm - S3
- 2.8mm / 4.5mm - S5
- 2.8mm / 4.5mm - S7



**AR INSTRUMED**

Dental Implant  
Zahnimplantat  
Implante dental  
Implant dentaire

## Conical Drills Universal Implant Kit



108-007

Conical Drills  
Universal Implant  
Kit

### 8pcs Trepbine Drills

- Int: 2.0mm, Ext: 3.0mm
- Int: 3.0mm, Ext: 4.0mm
- Int: 4.0mm, Ext: 5.0mm
- Int: 5.0mm, Ext: 6.0mm
- Int: 6.0mm, Ext: 7.0mm
- Int: 7.0mm, Ext: 8.0mm
- Int: 8.0mm, Ext: 9.0mm
- Int: 9.0mm, Ext: 10.0mm



108-008

Conical Drills  
Universal Implant  
Kit

### 8pcs Titanium Coated Drills

- Dia: 2.0mm, L: 18.5mm
- Dia: 2.5mm, L: 18.5mm
- Dia: 2.8mm, L: 18.5mm
- Dia: 3.2mm, L: 18.5mm
- Dia: 3.65mm, L: 18.5mm
- Dia: 4.5mm, L: 18.5mm
- Dia: 4.8mm, L: 18.5mm
- Dia: 5.2mm, L: 18.5mm



# AR INSTRUMED

Dental Implant  
Zahnimplantat  
Implante dental  
Implant dentaire

## Conical Drills Universal Implant Kit



108-009

Conical Drills  
Universal Implant  
Kit

### The Kit Includes:

- Compression Screw 2.8 / 4.0mm
- Compression Screw 4.0mm
- Driver 2.5mm
- Compression Screw 2.0 / 2.8mm
- Compression Screw 2.8 / 3.4mm
- Tissue Punch 4.5mm
- Tissue Punch 5.8mm
- Pilot Drill 2.0mm
- Drill Extender
- Ratchet Driver 2.5mm
- Hand Driver 2.5mm
- Torque Wrench (10-40Ncm)
- Depth Gauge

## Implant Trimmer Kit



108-010

Implant Trimmer Kit  
Kit

### The Kit Includes:

- |                 |              |             |           |               |
|-----------------|--------------|-------------|-----------|---------------|
| Trephine Drills | Tissue Punch | Guide Punch | Saw Disks | Sterilization |
| - 3.0mm         | - 3.5mm      | - 4.0mm     | - 7.0mm   | Cassette/ Box |
| - 4.0mm         | - 4.0mm      | - 5.0mm     | - 10.0mm  |               |
| - 5.0mm         | - 4.5mm      |             |           |               |
| - 6.0mm         | - 5.0mm      |             |           |               |
| - 7.0mm         | - 5.5mm      |             |           |               |



# AR INSTRUMED

Dental Implant  
Zahnimplantat  
Implante dental  
Implant dentaire

## Mini Multi Hex Drivers Kit



108-011

Mini Multi Hex  
Drivers Kit

A Unique designed Mini Hex Drivers Kit, includes the widely used Hex Drivers during an implant Surgery with a Ratchet Extender & an Infinite Torque Ratchet Wrench.

### The Kit Includes:

- Ratchet Hex Drivers 1.25mm Long & Short
- Ratchet Hex Drivers 2.42mm Long & Short
- Ratchet Extender 4.0mm
- Hand Hex Driver 1.25mm Long
- Infinite Torque Ratchet Wrench 4.0mm Square Driver
- Smart Sterilization / Organizer Kit

## Internal Irrigation Drills Kit



108-012

Internal Irrigation  
Drill Kit

The use of internal irrigation drilling systems is promoted on the basis of improved delivery of coolant to the bone/drill interface, which should subsequently improve local debridement and cooling, and therefore reduce thermal insults to the bone thus limiting frictional heat and reducing risk of osteonecrotic changes caused due to overheating.

### The Kit Includes:

#### 8pcs Internal Irrigation Drills

- Dia: 2.0mm, L: 18.5mm
- Dia: 2.4mm, L: 18.5mm
- Dia: 2.8mm, L: 18.5mm
- Dia: 3.2mm, L: 18.5mm
- Dia: 3.6mm, L: 18.5mm
- Dia: 4.2mm, L: 18.5mm
- Dia: 4.6mm, L: 18.5mm
- Dia: 5.0mm, L: 18.5mm



**AR INSTRUMED**

Dental Implant  
Zahnimplantat  
Implante dental  
Implant dentaire

### Trephine Drills Kit



108-013

Trephine Drills  
Kit

The trephines are used to make osseous takings with the aim of an autogenous transplant. The sharpened blades are to increase the cutting efficiency to realize a non traumatic treatment.

The depth markings make it easier to control the Drilling Depths efficiently.

#### 8pcs Trephine Dimensions

Internal Ø	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm	7.0mm	8.0mm	9.0mm
External Ø	3.0mm	4.0mm	5.0mm	6.0mm	7.0mm	8.0mm	9.0mm	10.0mm

### Implant Tissue Punch Kit



108-014

Implant Tissue  
Punch Kit

A tissue punch is a simple and easy way to gain access to the underlying bone.

It generally is chosen in Dental practice when there is:

- 1) Excellent ridge width and there will be no bone needed i.e. mandibular posterior,
- 2) Excellent keratinized gingiva,
- 3) Good access to place the punch.

#### 6pcs Latch Type Tissue Punch

- Dia: 3.5mm
- Dia: 4.0mm
- Dia: 4.5mm
- Dia: 5.0mm
- Dia: 5.5mm
- Dia: 6.0mm



**AR INSTRUMED**

Dental Implant  
Zahnimplantat  
Implante dental  
Implant dentaire

### Conical Drills Kit



**25pcs Conical Drills with Stoppers**

- 1.8mm / 2.2mm - S1
- 1.8mm / 2.2mm - S2
- 1.8mm / 2.2mm - S3
- 1.8mm / 2.2mm - S5
- 1.8mm / 2.2mm - S7
- 2.2mm / 2.7mm - S1
- 2.2mm / 2.7mm - S2
- 2.2mm / 2.7mm - S3
- 2.2mm / 2.7mm - S5
- 2.2mm / 2.7mm - S7
- 2.3mm / 3.3mm - S1
- 2.3mm / 3.3mm - S2
- 2.3mm / 3.3mm - S3
- 2.3mm / 3.3mm - S5
- 2.3mm / 3.3mm - S7
- 2.8mm / 3.7mm - S1
- 2.8mm / 3.7mm - S2
- 2.8mm / 3.7mm - S3
- 2.8mm / 3.7mm - S5
- 2.8mm / 3.7mm - S7
- 2.8mm / 4.5mm - S1
- 2.8mm / 4.5mm - S2
- 2.8mm / 4.5mm - S3
- 2.8mm / 4.5mm - S5
- 2.8mm / 4.5mm - S7

### Tri-Spade Implant Drills (Polished)



104-016	104-017	104-018	104-019	104-020	104-021	104-022	104-023
2mm	2.8mm	3.2mm	3.65mm	4.2mm	4.5mm	4.8mm	5.2mm

High Quality Stainless Steel. Highly tempered, Black Titanium NiTride coated for much Longer Life. Drills remain Sharper for Prolonged Period of Time as compared to the Non-Coated Drills.

#### Drill Dimensions:

- Dia: 2.0mm, L: 18.5mm
- Dia: 2.8mm, L: 18.5mm
- Dia: 3.2mm, L: 18.5mm
- Dia: 3.65mm, L: 18.5mm
- Dia: 4.2mm, L: 18.5mm
- Dia: 4.5mm, L: 18.5mm
- Dia: 4.8mm, L: 18.5mm
- Dia: 5.2mm, L: 18.5mm

The tri-spade drill design, which features three cutting edges, is much more stable in the hands of the clinician.

It reduces crestal chatter upon entry into the bone site (a stable drilling situation) resulting in a more perfectly prepared final hole for placement of the dental implant.



**AR INSTRUMED**

Dental Implant  
Zahnimplantat  
Implante dental  
Implant dentaire

### Titanium NiTride Coated Drills



104-024	104-025	104-026	104-027	104-028	104-029	104-030	104-31
2mm	2.8mm	3.2mm	3.65mm	4.2mm	4.5mm	4.8mm	5.2mm

High Quality Stainless Steel, Highly tempered, Black Titanium NiTride coated for much Longer Life,  
Drills remain Sharper for Prolonged Period of Time as compared to the Non-Coated Drills.

#### Drill Dimensions:

- Dia: 2.0mm, L: 18.5mm
- Dia: 2.8mm, L: 18.5mm
- Dia: 3.2mm, L: 18.5mm
- Dia: 3.65mm, L: 18.5mm
- Dia: 4.2mm, L: 18.5mm
- Dia: 4.5mm, L: 18.5mm
- Dia: 4.8mm, L: 18.5mm
- Dia: 5.2mm, L: 18.5mm

### Tri-Spade Implant Drills with Titanium NiTride Coating



104-032	104-033	104-034	104-035	104-036	104-037	104-038	104-39
2mm	2.8mm	3.2mm	3.65mm	4.2mm	4.5mm	4.8mm	5.2mm


High Quality Stainless Steel, Highly tempered, Gold Titanium Nitride coated for much Longer Life,  
Drills remain Sharper for Prolonged Period of Time as compared to the Non-Coated Drills.

#### Drill Dimensions:

- Dia: 2.0mm, L: 18.5mm
- Dia: 2.8mm, L: 18.5mm
- Dia: 3.2mm, L: 18.5mm
- Dia: 3.65mm, L: 18.5mm
- Dia: 4.2mm, L: 18.5mm
- Dia: 4.5mm, L: 18.5mm
- Dia: 4.8mm, L: 18.5mm
- Dia: 5.2mm, L: 18.5mm



A vertical ruler scale is positioned on the right side of the image, ranging from 0 to 20. The scale marks are clearly visible, providing a reference for the size of the instruments shown in the photographs.





A vertical ruler scale is positioned on the right side of the image, showing measurements from 0 to 20 in increments of 1. The scale is oriented vertically along the right edge of the frame.



# AR INSTRUMED DEUTSCHLAND GMBH

+49 30 91201490

[www.ar-instrumed.de](http://www.ar-instrumed.de)

AR INSTRUMED Deutschland GmbH  
Buttmann Str. 7A, 13357 Berlin Germany  
Tel: +49 30 91201490-1 Fax: +49 30 91201492  
[www.ar-instrumed.de](http://www.ar-instrumed.de) - [info@ar-instrumed.de](mailto:info@ar-instrumed.de)

Distributed By